


IMPLANTOLOGY

 Outlink²




sweden & martina

IMPLANTOLOGY TABLE OF CONTENTS



Implant system

06



The range

18



Prosthetic components

46



Surgical instruments

20



Material composition


86



Guide to the choice
of prosthetic solutions

44

Sweden & Martina develops and manufactures implant systems that offer both excellent clinical functionality and perfect aesthetic results. The surfaces have been designed to obtain the best relationships between the surface roughness of the titanium and speed of bone healing. The surgical instruments are functional, easy to use and ergonomic. Training courses, continuous refresher courses and extensive assistance distinguish the service and reliability that have made Sweden & Martina a leader in the Italian implant market.




Multifunctional mounter:

The particular conformation of the Outlink² mounter allows it not only to act as a carrier for the transport and positioning of the implant in its site, but also as a transfer for taking the impression and as post.



The Outlink² implant system is suitable for:

- standard operating procedures involving the double or single surgical phase;
- immediate loading;
- post-extraction situations;
- guided regeneration protocols;
- Switching Platform.



Outlink² Shorty:

The Outlink² implant range also includes a line of short fixtures, ideal in the case of small vertical bone dimension.



Wide range of fixtures:


The Outlink² implant system has a complete range of fixtures for diameters and heights, thus offering extremely versatile implant-prosthetic solutions. The Outlink² implants are available in diameters 3.30, 3.75, 4.10 (with the two different connections, standard and SP) and 5.00, and in all heights from h 5 mm to h 18 mm.





Cylindrical shape:

The implant body is cylindrical with a conical apical conformation for greater ease of insertion, and its external hexagon connection allows it to be used in different clinical situations.



Switching Platform:

The geometrical characteristics of Outlink² connections enable them to be applied in Switching Platform protocols.



Completeness of prosthetic solution:

Ideal for both cemented and screwed prostheses, the external hexagon connection makes it particularly suitable in situation of severe disparallelism.




Scientific support:

The use of the Outlink² system is backed up by numerous publications which document its extreme versatility with over 10 years of clinical success.


Outlink² ø 3.30 implant

The Outlink² ø 3.30 implant has a platform with 2.40 mm external hexagon, 1.0 mm high and threading M 1.8, allowing any type of prosthetic restoration to be produced satisfactorily.

Due to their small diameter, ø 3.30 implants are ideal for implant-prosthetic rehabilitation when there is limited space between adjacent teeth, as in the case of single crowns in upper lateral incisor positions and in lower intraforaminal situations.*





Stability with respect to disto-mesial and antero-posterior stress guaranteed by the external hexagon height of 1.00 mm. Resistant and particularly sturdy section despite its small diameter, thanks to the external hexagon connection.



The thread of the implants has a pitch of 0.6 mm which facilitates screwing progress and limits bone trauma after application of the load.



Tapered apex with large discharge notches that give the implant excellent self-tapping properties; the fully threaded apex section considerably simplifies its insertion.



* They can also be used for the rehabilitation of single crowns at premolar level. In distal sectors they must be used exclusively for the rehabilitation of multiple fixed structures. They are also very useful in the case of total edentulism on thin mandibular crests where it is preferred not to carry out regeneration. In this case it is recommended to use at least 4 fixed implants with a bar.



Outlink² ø 3.75 and ø 4.10 implants

The Outlink² implant with prosthetic platform ø 4.10, with 2.70 mm standard hexagon 0.70 mm high and threading of M 2.0, is available both with a 4.10 mm neck and 3.75 mm spire and with a 4.10 mm neck and 4.10 mm spire.

Using the same platform (4.10 mm) it is possible to choose between two different spire diameters, 3.75 mm and 4.10 mm, depending on the available bone thickness.




The connection platform of Outlink² implants has an external hexagon that today is generally recognised as standard at world level.

The external connection makes them particularly suitable for operations in the case of multiple edentulism with severe disparallelism, as it considerably facilitates the phase of taking the impression and the subsequent insertion and removal of prostheses.

Outlink² ø 4.10 SP implant (Switching Platform)

The Outlink² ø 4.10 SP implant has a 4.10 mm prosthetic platform, a 2.40 mm hexagon 1.0 mm high with a threading of M 1.8, the same as those of the ø 3.30 implant. This characteristic allows the use of prosthetic components with diameter 3.30 mm, optimally performing the Switching Platform technique which takes advantage of the horizontal component of the biological width, thus minimising the loss of crestal bone.




The Switching Platform is a prosthetic rehabilitation technique that requires the use of posts with a smaller diameter than the implant platform in order to improve the biomechanical distribution of the prosthetic load, but especially to distance the prosthetic connection from the cervical bone.

The portion of the connection platform not occupied by the prosthesis creates a supporting base for the connective tissue, thus stabilising the collagen fibres and in this way minimising bone reabsorption.



The Switching Platform technique is possible with ø 4.10 mm SP implants, using ø 3.30 mm prosthetic components on these implants.




FEM analysis of Outlink² ø 4.10 mm implants with a ø 3.30 mm post according to the Switching Platform protocol.




Outlink² ø 5.00 implant

The Outlink² ø 5.00 presents a prosthetic platform with diameter 5.00 mm with 2.70 mm external hexagon, 0.70 mm high and threading M.20, the same as those of the standard 4.10 mm platform, which guarantees high precision and versatility. Outlink² ø 5.00 implants allow the application of the Switching Platform technique using ø 4.10 mm prosthetic components.



The 5.00 mm diameter of this implant makes it ideal for implant-prosthetic rehabilitation on thick bone crests. The external connection and the wide diameter of the spires give this implant extraordinary sturdiness and stability.



FEM analysis of Outlink² ø 5.00 mm implants with a ø 4.10 mm post according to the Switching Platform protocol.



The Switching Platform technique is possible with ø 5.00 mm implants, using ø 4.10 mm prosthetic components on these implants. This improves the preservation of the crestal bone.

Multifunctional mounter

The Outlink² implant has the mounter already assembled in the PMMA vial. As well as the traditional carrier function for the transport and positioning of the implant, the particular conformation of the Outlink² mounter, made in Gr.5 titanium, also allows it to be used as a transfer when taking the impression and as a post during prosthetic rehabilitation.



The retentive tabs can be cut easily, so that the mounter adapts to the morphology of the element to be prosthodontized and thus become a practical post.

The conical profile of the mounter facilitates insertion and removal of the crowns or of bridges, in the case of multiple structures.


The mounter is supplied already preassembled with the implant. The connecting screw is also available separately as a spare

The golden colour of the mounter/post guarantees maximum results as regards the aesthetic appearance of reconstructions.



..... The thickness of the mounter is such as to allow it to be reduced in height if necessary, or milled, and to create coulisses in the walls for repositioning the prosthesis.

The mounters with all diameters have two repositioning faces to guarantee a good non-rotational aspect while taking the impression.



Practicality of the surgical procedure: the view of the mounter from above shows the conformation of the upper part, with an internal octagon, which allows it to be easily lifted by the driver and put into position.



Outlink² Shorty implants

Shorty Outlink² fixtures with height 5.0, 7.0 and 8.5 mm are available in the program; they can be used, according to the most recent clinical protocols, in all cases where there is small vertical bone dimension.



WARNING: Never use these implants for rehabilitating single crowns, but only as support posts combined with longer fixtures for multiple rehabilitations. It is also recommended to always use, whenever possible, implants with the largest diameter possible depending on the thickness of the crest.



In case of very short implants (5 mm and 7 mm), the apical tapering was redesigned to improve primary stability even further. Installing a prosthesis with the Switching Platform technique is recommended for these implants in order to preserve the already reduced vertical dimension of the crest as much as possible. This choice is necessary in Outlink² Shorty implants with a diameter of 4.10 mm because they have a 2.4 mm hexagon instead of the 2.7 mm standard hexagon (4.10SP platform).



Key to the implant codes

The implant codes are so-called "talking" codes, i.e. they allow easy identification of the piece. Below is a table showing how the talking codes work using Code E2-ZT-410-SP-115 as an example:


| Type of implant E2- E2- | Surface ZT- ZT- | Diameter 410 410 | Connection SP SP | Length 115 115 |
|----------------------------------|---|--|---|---|
| E2: Outlink ² implant | 3S: Trisurface Surface ZT: ZirTi surface | 330: 3.30 mm 375: 4.10 mm 410 - 410SP: 4.10 mm 500: 5.00 mm | SP: Switching Platform(e.g. 2.40 mm) | 050: 5 mm 070: 7 mm 085: 8.5 mm 100: 10 mm 115: 11.5 mm 130: 13 mm 150: 15 mm |
| | | <i>is the size of the diameter of the implant connection</i> | | <i>refers to the length of the implant</i> |

Table of colour codes

A colour code system has been defined in the Outlink² implant system for identifying the endosseous diameter of the implant (see table on page 44-45).

The colour code identifies:


- the transfers for taking an impression and the laboratory analogs;
- the final drills;
- the sequence on the surgical tray.

| | ø 3.30 | ø 3.75 | ø 4.10 | ø 4.10SP | ø 5.00 |
|-------------------------|---|---|---|---|---|
| Colour code on the pack |  |  |  |  |  |


Surfaces

These studies have shown that the closer the roughness is to the size of the osteoblasts, the more influence it has on cell behaviour, causing the platelet activity to increase with respect to a smooth surface, thus accelerating the repair and osseointegration processes. The roughness is able to orient the cell layout, to influence their metabolism and proliferation, to differentiate osteoblasts and to modulate the production of extra-cellular matrix. These studies have led to the current development of the Outlink² implant surfaces: ZirTi (Zirconium Sand Blasted Acid Etched Titanium) and TriSurface.


ZirTi surface



This is a surface in which the roughness is obtained with subtraction techniques by sand-blasting with zirconium oxide and acid-etching with mineral acids.




The roughness assumed by the surface of the implant body is an ideal situation for promoting osteoblastic proliferation and differentiation, as well as the formation and growth of bone tissue.




Histology showing vital mineralised bone in intimate contact with the ZirTi surface of an Outlink² implant. Osteons and areas of bone marrow can be seen.

Image and caption by kind permission of Dr Daniele Botticelli.




In both types of surface the collar is machined for 0.75 mm of the height.


TriSurface Surface



The implant with a TriSurface surface has a sandblasted coronal portion, in order to obtain an intermediate level of roughness that allows better control of any bacterial infections before they can degenerate into peri-implantitis.



The middle apical portion of the body of the implant is coated with HRPS (High Roughness Plasma Spray) and has the maximum level of roughness that can be obtained, thus guaranteeing excellent primary stability even when the bone is only slightly mineralised, and significantly increasing the bone-implant contact surface.



Histological image of the bone growth around an Outlink² implant with TriSurface surface.

*Il Circolo
Rivista Periodica di Odontostomatologia,
1: 13-20, 2004*

Cold plasma surface decontamination

The better the processes of passivation, cleaning and decontamination of an implant surface, the greater the presence of pure titanium able to come in contact with the bone. This proportionally increases the possibilities of osseointegration.


At the end of the surface treatments, the implants are subjected to a careful cleaning and decontamination process by means of cold plasma triggered in Argon after first being cleaned of the main processing residue with numerous washing cycles in specific solvents.

What is decontamination? It is the total removal of dirt, particle residue and Bioburden from the surface of the implants, carried out before sterilisation.

During the Argon treatment, the gas atoms are partially ionised, they acquire energy and "bombard" the surface of the fixture violently. This kind of "atomic sand-blasting" removes organic contaminants without leaving any traces or additional residuals.



Implant before the decontamination treatment.



Working plasma reactor during surface decontamination of the implants.



Implant after the decontamination treatment.

As known, Argon is an inert gas that does not react with the titanium surfaces. The condition of surface decontamination is controlled regularly with randomised analyses of Bioburden residuals and a SEM visual examination on all the batches produced. This process activates the ionisation of the atoms on the surface of the titanium oxide which in turn increases the wettability of the fixture.



Implant packaging

The implants are packaged in PMMA vials in which they are held/housed in special titanium baskets that protect the surface of the fixture against possible recontaminations due to contact. All the materials comprising the packaging have been suitably tested to verify their suitability to sterilisation, preservation and medical use.



The blister packs containing the implants are packed in cardboard boxes which also contain the labels for the patient records reporting the details for traceability of the product (code and batch number).



The vials are contained in a special PETG blister pack sealed by a Tyvek cover, which guarantees sterility of the product for 5 years.



The mounter, assembled on every single implant, is visible, ready to be engaged by the special drivers.

The surgical cover screw, supplied with every single implant, is housed in a special compartment in the top of a blue cap made of LDPE which closes the vial. In turn, a small transparent lid in PMMA closes the blue cap.

Sterilisation











Sterilisation is the total elimination of the residual microbial load present on the implant after the decontamination and packing process, it is carried out with the use of beta rays. The sterilisation procedures are carried out in accordance with the UNI EN ISO 13485 and UNI EN ISO 9001 quality standards.

A beta ray sterilisation process was chosen because it has a variety of different advantages:

- the process occurs in a completely automatic way with computerised control of all the phases;
- the process is quick, reliable and extremely easy to repeat with safety and precision;
- the process is extremely eco-friendly, does not require the presence of radioactive sources and does not lead to the formation of toxic or radioactive products;
- beta rays are minimally invasive with regards to packaging due to the speed of the treatment. This guarantees preservation of the product's sterility over time (certified duration of 5 years).

THE RANGE









Outlink² implants

| implant diameter | ø 3.30 mm | ø 3.75 mm | ø 4.10 mm | ø 4.10SP mm | ø 5.00 mm |
|-------------------------|---|---|---|---|---|
| mounter* |  E-MOU2-330 |  E-MOU2-410 |  E-MOU2-410 |  E-MOU2-330 |  E-MOU2-500 |
| connecting screws |  VM-180 |  VM-200 |  VM-200 |  VM-180 |  VM-200 |
| surgical cover screws** |  E-VT-330 |  E-VT-410 |  E-VT-410 |  E-VT-330 |  E-VT-500 |

* The mounters are sold preassembled with the implants. Both the mounters and the connecting screws (VM-180 and VM-200) are available on sale as individual spare parts. If the mounter is used as a post, the torque for tightening the screws is 20-25 Ncm.







** Each implant is sold with its own surgical cover screw. The surgical screws are also available on sale individually in a sterile pack and must be tightened to 10 Ncm.

Shorty implants

| implant diameter | ø 3.30 mm | ø 3.75 mm | ø 4.10 mm | ø 4.10SP mm | ø 5.00 mm |
|---------------------|-----------|--|--|--|--|
| 5.00 | - | - | - |  5.00 |  5.00 |
| ZirTi TriSurface | | | | E2-ZT-410SP-050 | E2-ZT-500-050 |
| 7.00 | - | - | - |  7.00 |  7.00 |
| ZirTi TriSurface | | | | E2-ZT-410SP-070 | E2-ZT-500-070 |
| 8.50 | - |  8.50 |  8.50 |  8.50 |  8.50 |
| ZirTi TriSurface | | E2-ZT-375-085 E2-3S-375-085 | E2-ZT-410-085 E2-3S-410-085 | E2-ZT-410SP-085 E2-3S-410SP-085 | E2-ZT-500-085 E2-3S-500-085 |




Standard implants

| implant diameter | ø 3.30 mm | ø 3.75 mm | ø 4.10 mm | ø 4.10SP mm | ø 5.00 mm |
|---------------------|---|---|---|---|---|
| 10.00 |  |  |  |  |  |
| ZirTi TriSurface | E2-ZT-330-100 E2-3S-330-100 | E2-ZT-375-100 E2-3S-375-100 | E2-ZT-410-100 E2-3S-410-100 | E2-ZT-410SP-100 E2-3S-410SP-100 | E2-ZT-500-100 E2-3S-500-100 |
| 11.50 |  |  |  |  |  |
| ZirTi TriSurface | E2-ZT-330-115 E2-3S-330-115 | E2-ZT-375-115 E2-3S-375-115 | E2-ZT-410-115 E2-3S-410-115 | E2-ZT-410SP-115 E2-3S-410SP-115 | E2-ZT-500-115 E2-3S-500-115 |
| 13.00 |  |  |  |  |  |
| ZirTi TriSurface | E2-ZT-330-130 E2-3S-330-130 | E2-ZT-375-130 E2-3S-375-130 | E2-ZT-410-130 E2-3S-410-130 | E2-ZT-410SP-130 E2-3S-410SP-130 | E2-ZT-500-130 E2-3S-500-130 |
| 15.00 |  |  |  |  | - |
| ZirTi TriSurface | E2-ZT-330-150 E2-3S-330-150 | E2-ZT-375-150 E2-3S-375-150 | E2-ZT-410-150 E2-3S-410-150 | E2-ZT-410SP-150 E2-3S-410SP-150 | - |
| 18.00 | - | - |  | - | - |
| ZirTi TriSurface | - | - | E2-ZT-410-180 | - | - |

Surgical kit

The Outlink² surgical kit has been designed and made to offer ease of use and immediate placing in the sequence of instruments*. The instruments, all made of stainless steel, have their codes screen-printed on the tray to allow the user to identify each instrument more easily and to put it back after the cleansing and cleaning phases, with the aid of a colour code system that traces the suitable surgical procedures for the various implant diameters. The Outlink² surgical kit is also supplied with the templates for the graphic representation of the implant measurements to allow choosing the most suitable implant diameters and lengths by means of radiographic or tomographic analyses






The kit consists of a practical box in Radel with a surgical tray inside that is set-up to hold the instruments according to a guided procedure. The sequences of use of the instruments are indicated by coloured marks.



A practical ratchet is also included that acts as a dynamometric key for checking the closing torque of the prosthetic screws and as a surgical key for inserting the implants. The ratchet has a very small head, making it easy to use even in distal sectors.

*Please note: to guarantee maximum duration of the surgical and prosthetic instruments, it is advisable to follow the recommended cleansing and sterilisation procedures.



| code | description |
|---|---|
|  | Surgical kit complete with all the instruments necessary for Outlink² implants |
| ZOUTLINK2* | |
|  | Radel instrument tray for Outlink² implants |
| OUT-KIT* | |
|  | Kit with 5 spare silicon supports for surgical trays, for drills or instruments with right-angle shanks |
| GROMMET-3 | |
|  | Kit with 5 spare silicon supports for surgical trays, for instruments fitted with connection hexagon |
| GROMMET-4 | |
|  | Kit with 5 spare silicon supports for surgical trays, for digital or handheld instruments |
| GROMMET-5 | |

* The abbreviations ZOUTLINK2* and OUT-KIT* are followed by a letter and number that indicate the revision of the kit. The contents of the Surgical Kit can be updated and varied according to the most effective and innovative surgical techniques.

OneBox²

The OneBox² surgical kit was created to meet the needs of surgeons who carry out a large number of implant operations and who therefore want to have a compact kit equipped essentially with all that is needed only for the surgical phase.





The OneBox² is a compact kit that is easy to carry, containing the surgical instruments strictly necessary for inserting Outlink2 implants.



The bone taps are present in the kit in the version with right-angle attachment.

This kit contains space for an additional series of alternative drills, which the surgeon may wish to choose and use depending on his experience.



| code | description |
|--|---|
|  ZEONEBOX* | OneBox ² surgical kit |
|  EONEBOX-KIT* | Instrument tray for OneBox ² |
|  GROMMET-CA-1 | Kit with 5 spare silicon supports for surgical trays, for drills or instruments with right-angle shanks |
|  GROMMET-CA-2 | Kit with 5 spare silicon supports for surgical trays, for instruments fitted with connection hexagon |

* The words ZEONEBOX* and EONEBOX-KIT* are followed by a letter and a number that indicate the revision of the kit. The contents of the kit can be updated and varied according to the most effective and innovative surgical techniques

WARNING: OneBox² does not contain drill depth stops or prosthetic drivers, but it contains all the drivers in the one-piece digital version and the right-angle version, which are much more practical during surgical procedures.

Screw Kit

The Sweden & Martina Screw Kit is a handy set containing the necessary drivers* for the prosthetic phases after removal of the Outlink² transgingival healing screws.

It includes digital and right-angle drivers, as well as a dynamometric ratchet.

Small and easy to carry, it allows simple and immediate management of the prosthetic rehabilitation phase after surgery.





As well as digital and right-angle drivers, the Screw Kit includes a carrier for offset P.A.D. abutments, thus also favouring rapid full-arch prosthetic rehabilitations.



The compact dimensions of this new kit and the intuitive layout of the tray with the names of the instruments indicated in their place, allow the surgeon practical and easy management of the post-surgical rehabilitation phase, with notable optimisation of the time spent at the chair.

The Screw Kit contains all the drivers for the various prosthetic solutions: for standard posts, for abutments, for P.A.D. prostheses, for Locator abutments, for ball attachments and the respective retainer caps.



| code | description |
|--|---|
|  <p>ZSCREW*</p> | Complete Screw Kit |
|  <p>SCREW-TRAY*</p> | Instrument tray for Screw Kit |
|  <p>GROMMET-CA-1</p> | Kit with 5 spare silicon supports for surgical trays, for drills or instruments with right-angle shanks |
|  <p>GROMMET-CA-2</p> | Kit with 5 spare silicon supports for surgical trays, for instruments fitted with connection hexagon |

* The abbreviations ZSCREW* and SCREW-TRAY* are followed by a letter and a number that indicate the revision of the kit. The contents of the kit can be updated and varied according to the most effective and innovative surgical techniques.

Initial, intermediate and countersink drills

All Sweden & Martina drills are made of stainless steel with high resistance to corrosion and wear. The extreme accuracy of design and production allows use completely free from vibrations and oscillations.

Precision drill: very sharp and precise, it facilitates the opening of the implant site especially in the case of very hard cortical bone.

Intermediate drills: the shape with progressive diameter of the 2 intermediate drills allows gradual expansion of the bone when using a drill with a diameter much larger than the previous one, so as to reduce bone stress to a minimum.




Pilot drill: the different heights of the laser markings allow better visibility in the surgical field and greater ease of recognition of the different depths during use.

Countersink drill: ideal for preparing the coronal part of the site in the case of implants with a prosthetic platform wider than the diameter of the spires.



initial drill, intermediate drills and countersink drill

kit



ZOUTLINK2*
ZEONEBOX*

FS-230
precision drill

FG-200/280XS
intermediate drill
ø 2.00 - 2.80 mm


FG-330/425XS
intermediate drill
ø 3.30 - 4.25 mm

FC-XS
countersink drill

pilot drill

pilot drill stops

kit



ZOUTLINK2*
ZEONEBOX*

FPT*-200-LXS
pilot drill

STOP200-070**
7.00 mm pilot
drill stop

STOP200-085**
8.50 mm pilot
drill stop

STOP200-100**
10.00 mm pilot
drill stop


STOP200-115**
11.50 mm pilot
drill stop

STOP200-130**
13.00 mm pilot
drill stop

STOP200-150**
15.00 mm pilot
drill stop

* The word STOP is followed by a number indicating the revision of the accessory. The stops are included only in the ZOUTLINK2* surgical kit.

Lt: Total length of the working part, including the tip.
Lp: Length of the tip. This measurement must be calculated in addition to the length of the preparation hole.



Please note: The letters FPT are followed by a number (2, 3) indicating the length of the drill shank: 2 indicates a length of 12.5 mm, 3 indicates a length of 14 mm. All the STOP2 and STOP3* are functional to any of these batches.

Please note: The initial drills (FPT*) and the final drills (FFT*, shown on the following page) always make a hole longer than the implant that is to be inserted. The oversizing (Lp) is equal to the height of the tip of the drill that is being used.

Final drills and stops

Also made of stainless steel with high resistance to corrosion and wear, Outlink² final drills present a number of cutting edges proportional to the hole diameter, so as to allow a continuous and homogeneous cutting movement and greater instrument stability during operation. All this enables high-precision implant preparations to be obtained, with consequent ease in inserting the implant.































The stops are inserted from the tips of the drills, with great ease and speed.

The coloured rings make it easy to recognise the dedicated instruments for each implant diameter.

The surgical drills have a spiral designed to allow the operator better control of the forward movement and better centring during the drilling phase.



| implant diameter | ø 3.30 mm | ø 3.75 mm | ø 4.10 mm ø 4.10SP mm | ø 5.00 mm | kit |
|------------------|--|--|--|--|-------------------------|
| drills** |  0.81 FFT*-280-LXS 19.60 |  0.87 FFT*-300-LXS 19.60 |  0.95 FFT*-340-LXS 19.70 |  1.23 FFT*-425-LXS 20.00 | ZOUTLINK2* ZEONEBOX* |
| |  STOP*-280-070 |  STOP*-300-070 |  STOP*-340-070 |  STOP*-425-070 | ZOUTLINK2* |
| |  STOP*-280-085 |  STOP*-300-085 |  STOP*-340-085 |  STOP*-425-085 | ZOUTLINK2* |
| |  STOP*-280-100 |  STOP*-300-100 |  STOP*-340-100 |  STOP*-425-100 | ZOUTLINK2* |
| stop |  STOP*-280-115 |  STOP*-300-115 |  STOP*-340-115 |  STOP*-425-115 | ZOUTLINK2* |
| |  STOP*-280-130 |  STOP*-300-130 |  STOP*-340-130 |  STOP*-425-130 | ZOUTLINK2* |
| |  STOP*-280-150 |  STOP*-300-150 |  STOP*-340-150 |  STOP*-425-150 | ZOUTLINK2* |

***Please note:** The letters FFT are followed by a number (2, 3) indicating the length of the drill shank: 2 indicates a length of 12.5 mm, 3 indicates a length of 14 mm. All the STOP2 and STOP3* are functional to any of these batches.

****Please note:** The drills always make a hole that is longer than the implant to be inserted. The oversizing (Lp) is equal to the height of the tip of the drill that is being used. See image on page 27.

Bone taps

Outlink² implants are self-tapping implants with excellent cutting and insertion capabilities. However, the use of a bone tap is recommended in all cases where the type of bone requires it. The absence of tapping in cases where this is recommended may lead to problems later when inserting the implant. They are available both with right-angle shank and with a connector for dynamometric ratchet.

Standard bone tap: to be used with the dynamometric ratchet, gives great sensitivity during use.

Short bone tap: to be used with the dynamometric ratchet, it is very useful in distal sectors.

Right-angle bone tap: for better control of the inserting axis and more uniform preparation.






| implant diameter | ø 3.30 mm | ø 3.75 mm | ø 4.10 mm ø 4.10SP mm | ø 5.00 mm | kit |
|-----------------------|-------------|-------------|--------------------------|-------------|------------|
| standard bone taps | | | | | ZOUTLINK2* |
| | E-MS-330 | E-MS-375 | E-MS-410 | E-MS-500 | |
| right-angle bone taps | | | | | ZEONEBOX* |
| | E-MS-330-CA | E-MS-375-CA | E-MS-410-CA | E-MS-500-CA | |
| short bone taps* | | | | | - |
| | E-MSC-330 | E-MSC-375 | E-MSC-410 | E-MSC-500 | |

* Please note: short bone taps are not included in any surgical kit.

Osteotomes

A set of steel osteotomes is available, useful for maxillary sinus floor elevation via the crestal bone and bone expansion protocols. The laser-etched codes on the handles show the diameter and height of the corresponding implant, so as to facilitate recognition of the correct surgical sequence.



The use of osteotomes is ideal in bone expansion procedures.



| code | description |
|-------------|--|
| E-OS-020-PP | Outlink ² Osteotome ø 0.20, flat tip |
| E-OS-090-PP | Outlink ² Osteotome ø 0.90, flat tip |
| E-OS-160-PC | Outlink ² Osteotome ø 1.60, concave tip |
| E-OS-200-PC | Outlink ² Osteotome ø 2.00, concave tip |
| E-OS-240-PC | Outlink ² Osteotome ø 2.40, concave tip |

Drills for distal sectors

As an option, shorter drills are available that are very practical in distal sectors with limited oral opening. They come in a wide range of diameters and are also useful for preparations in extremely compact bone where, in the most coronal portion, you want to widen the preparation diameter by 0.10 mm with respect to the size of the standard drills to facilitate the insertion of the implants. On the other hand, in low-density bone they can be used to under-prepare the implant site so as to obtain optimum primary stability.



The drills have etched depth lines from 7 to 15 mm.






Cylindrical drills

ø 2.00 mm



FFT5-200-LXS*
Cylindrical drill ø 2.00 mm

ø 2.80 mm




FFT5-280-LXS*
Cylindrical drill ø 2.80 mm

ø 2.90 mm




FFT5-290-LXS*
Cylindrical drill ø 2.90 mm

ø 3.00 mm




FFT5-300-LXS*
Cylindrical drill ø 3.00 mm

ø 3.20 mm



FFT5-320-LXS*
Cylindrical drill ø 3.20 mm

ø 3.30 mm



FFT5-330-LXS*
Cylindrical drill ø 3.30 mm

ø 3.40 mm



FFT5-340-LXS*
Cylindrical drill ø 3.40 mm

ø 3.60 mm




FFT5-360-LXS*
Cylindrical drill ø 3.60 mm

ø 4.25 mm




FFT5-425-LXS*
Cylindrical drill ø 4.25 mm

ø 4.45 mm












FFT5-445-LXS*
Cylindrical drill ø 4.45 mm

* The drills for distal sectors are without irrigation and are not included in any surgical kit. They cannot be used with depth stops












Please note: The drills always make a hole that is longer than the implant to be inserted. The oversizing (L_p) is equal to the height of the tip of the drill that is being used.











Drivers and Screwdrivers

| code | description | kit |
|--|--|------------------------------------|
|  U-AW3-MOU-CA | Octagonal driver with right-angle attachment for Outlink ² implant mounter | ZOUTLINK2* |
|  U-AWV-MOUC | Short octagonal driver for Outlink ² implant mounter | ZOUTLINK2* |
|  AVV-CA-DG-EX | Hand knob for hand use of drivers, bone taps and drivers with right angle shank and with hexagonal connection for torque-control ratchet | ZEONEBOX* |
|  HSM-20-EX | Driver for connecting screws, with connector for dynamometric ratchet or digital connector, short | ZOUTLINK2* ZSCREW* |
|  HSML-20-EX | Driver for connecting screws, with connector for dynamometric ratchet or digital connector, long | ZOUTLINK2* ZSCREW* |
|  HSMXL-20-EX | Driver for connecting screws, with connector for dynamometric ratchet or digital connector, extra long | ZSCREW* |
|  HSMXS-20-DG | Driver for surgical connecting screws, digital, extra short | ZOUTLINK2* ZEONEBOX* ZSCREW* |
|  HSM-20-DG | Driver for connecting screws, digital, short | ZOUTLINK2* ZEONEBOX* ZSCREW* |
|  HSML-20-DG | Driver for connecting screws, digital, long | ZOUTLINK2* ZEONEBOX* ZSCREW* |












| code | description | kit |
|--|--|------------------------------------|
|  HSMXS-09-DG | Driver for surgical cover screws, digital, extra short | ZOUTLINK2* ZEONEBOX* |
|  HSM-09-DG | Driver for surgical cover screws, digital | ZOUTLINK2* ZEONEBOX* |
|  HSM-20-CA | Driver for connecting screws, with right-angle shank | ZOUTLINK2* ZEONEBOX* ZSCREW* |
|  HSM-09-CA | Driver for surgical cover screws, with right-angle shank | ZOUTLINK2* ZEONEBOX* |
|  AVV2-ABUT | Driver for standard abutments and for straight P.A.D. abutments | ZSCREW* |
|  BASCC-EX | Driver for ball attachments, with connector for dynamometric ratchet or digital connector | ZSCREW* |
|  8926-SW | Short driver in Gr. 5 titanium for screwing the Locator Abutments. The driver is compatible with the Outlink² system's dynamometric ratchet | ZSCREW* |
|  8927-SW | Long driver in Gr. 5 titanium for screwing the Locator Abutments. The driver is compatible with the Outlink² system's dynamometric ratchet | ZSCREW* |
|  PAD-CAR | Carrier for transferring angled P.A.D. abutments into the oral cavity, sterilisable and reusable. It must be fixed to the abutments with the screw PAD-VTRAL-140 | ZSCREW* |

Accessories

| code | description | kit |
|--|--|------------------------------------|
|  B-AVV-CA3 | Mechanical adapter with right-angle shank for instruments with hexagonal connector | ZOUTLINK2* ZEONEBOX* ZSCREW* |
|  BPM-15 | Extension for bone taps, mounters, drivers and manual drivers, with hexagonal connector for dynamometric key | ZOUTLINK2* ZEONEBOX* |
|  AVV3-MAN-DG | Hand knob for bone taps, mounters, drivers and manual drivers | ZOUTLINK2* ZSCREW* |
|  ORING180-088 | Kit with 5 spare o-rings for all accessories with hexagonal connector for dynamometric key | - |
|  CRI5-KIT | Kit composed of a ratchet, which can be used in dynamometric or fixed mode, and accessories for quick torque adjustment and periodic maintenance (Allen wrench and lubricant). The ratchet has torque limits from 10 to 70 Ncm, with adjustment lines at 10-20-25-30-35-50-70 Ncm. (Supplied with the Outlink ² surgical kit) | ZOUTLINK2* ZEONEBOX* ZSCREW* |
|  PP-2/28 | Parallelism pin Ø 2.00 and 2.80 mm | ZOUTLINK2* ZEONEBOX* |
|  PROF3 | Depth gauge | ZOUTLINK2* |
|  CMD | Mounter stop key | ZOUTLINK2* ZEONEBOX* |
|  E2-CM | Mounter stop key | ZOUTLINK2* ZEONEBOX* |
|  PROF-CAL2 | Extension for surgical drills | ZOUTLINK2* ZEONEBOX* |



| code | description | kit |
|--|---|-------------------------|
|  E-MOU-330 | Short mounter Ø 3.30 | ZOUTLINK2* |
|  E-MOU-410 | Short mounter Ø 4.10 | ZOUTLINK2* |
|  E-MOUL-330 | Long mounter Ø 3.30 | ZOUTLINK2* |
|  E-MOUL-410 | Long mounter Ø 4.10 | ZOUTLINK2* |
|  E-PAD-PS410-L | Bone profiler for levelling irregular bone crest, with wide flaring | - |
|  E-PAD-PS410-S | Bone profiler for levelling irregular bone crest, with narrow flaring | - |
|  E-L100 | X-ray template for Outlink² implants (real dimensions) | ZOUTLINK2* ZEONEBOX* |
|  E-L120 | X-ray template for Outlink² implants (dimensions increased by 20%) | ZOUTLINK2* ZEONEBOX* |
|  E-L130 | X-ray template for Outlink² implants (dimensions increased by 30%) | ZOUTLINK2* ZEONEBOX* |

Shorty drilling kit

The implant site of the specific Shorty implants with lengths of 5.0, 7.0 and 8.5 mm can be prepared with the drills contained in the standard surgical kits; however, the drills in these kits require an over-preparation connected to the measurement of the drill tip. On the other hand, the choice of a short implant is generally connected to the lack of available bone height in the implant site, so it would be preferable not to have to engage a working thickness with the tip of the drill, but instead to lodge a longer implant. For this reason, a drilling kit was created to enable preparing the sites for Shorty implants with a height of 5.7, and 8.5 mm and very short tip where a portion of over-preparation is not considered necessary. The Shorty drills allow dedicating the entire bone available to lodging the implant, without any waste. They also have the advantage of being shorter than the standard drills (24.85 mm instead of 35 mm). This important feature makes it possible to use these instruments even in case of difficult to reach distal sectors or small oral openings.




The kits contain two short parallelism pins, which are very handy in distal sectors.



All the instruments in the Shorty Drilling Kit are also available individually as spare parts.

The colour codes of the implant diameters facilitate the choice of the instruments.




| code | description |
|---|---|
|  | Drilling kit for short implants |
|  | Empty tray in Radel |
|  | Kit with 5 spare silicon supports for surgical trays, for drills or instruments with right-angle shanks |

* The abbreviations ZSHORTY* and SHORTY-KIT* are followed by a letter and a number that indicate the revision of the kit. The contents of the kit can be updated and varied according to the most effective and innovative surgical techniques.

Instruments contained in the Shorty Drilling Kit

Shorty pilot drill



FPS-200
Shorty pilot drill

Shorty pilot drill stops



STOPS-200-050
5.00 mm stop
for Shorty pilot drill



STOPS-200-060
6.00 mm stop
for Shorty pilot drill



















STOPS-200-070
7.00 mm stop
for Shorty pilot drill


| code | description | diameter |
|--|---|----------------|
|  FGS-200/300 | Short guide drill | 2.00/2.50/3.00 |
|  FGS-340/425 | Short guide drill | 3.40/3.80/4.25 |
|  FGS-425/540* | Short guide drill | 4.25/4.80/5.40 |
|  PPS-2/3 | Parallelim pin for short drills with depth lines at 5.00, 6.00 and 7.00 | 2.00/3.00 |

* Intermediate drill dedicated to other Sweden & Martina implant systems.



| | ø 3.75 mm | ø 4.10 mm ø 4.10SP mm | ø 5.00 mm | ø 6.00 mm |
|--------|--|--|--|---|
| drills |  24.85 FFS-300 |  24.85 FFS-340 |  24.85 FFS-425 |  24.85 FFS-540* |
| |  STOPS-300-050 |  STOPS-340-050 |  STOPS-425-050 |  STOPS-540-050* |
| stop |  STOPS-300-060 |  STOPS-340-060 |  STOPS-425-060 |  STOPS-540-060* |
| |  STOPS-300-070 |  STOPS-340-070 |  STOPS-425-070 |  STOPS-540-070* |

* Drills and stops with diameter 5.40 mm for other Sweden & Martina implant systems are also available in the Drilling Kit.



Please note: remember that the drills in the Drilling Kit do not over-prepare the surgical site. The working lengths include the portion related to the conical tip of the drill.

GUIDE TO THE CHOICE OF PROSTHETIC SOLUTIONS

implant diameter

Ø 3.30 mm

Ø 3.75 mm



label with colour code on the pack



final drills



FFT*-280-LXS



FFT*-300-LXS

connecting screw



VM-180



VM-200

open tray transfer



E-TRA-330-RIT



E-TRAR-330-RIT



E-TRAR-410-RIT

screw for open tray transfer



VTRA-180



VTRA-200

closed tray transfer




E-TRARS-330



E-TRARS-410

screw for closed tray transfer



E-VTRAS-180



E-VTRAS-200

analog



E-ANA-330



E-ANA-410



Ø 4.10 mm

Ø 4.10SP mm

Ø 5.00 mm



FFT*-340-LXS



FFT*-340-LXS



FFT*-425-LXS



VM-200



VM-180



VM-200



E-TRAR-410-RIT



E-TRAR-330-RIT



E-TRAR-410-RIT



E-TRAR-500-RIT



VTRA-200



VTRA-180



VTRA-200



E-TRARS-410




E-TRARS-330



E-TRARS-410



E-VTRAS-200



E-VTRAS-180



E-VTRAS-200



E-ANA-410



E-ANA-330



E-ANA-410



E-ANA-500

Transgingival healing caps

The transgingival healing caps in Gr. 5 titanium can be identified via a laser-marking that reports the diameter, emergence profile and height. In case of transgingival healing screws with straight emergence profile, the marking only reports the platform diameter and height. The transgingival healing caps must be tightened at a maximum torque of 10 Ncm, using the HSM series of drivers, the full details and codes of which can be found on page 36-37.

The laser marking indicates the diameter (in the example 33=3.30mm) and the transgingival height (in the example 5=5mm).



The laser marking indicates the connection diameter (in the example 50=5.00mm), the maximum crown size (in the example 60 = 6.00 mm) and the transgingival height (in the example 2 = 2.00 mm).

The transgingival caps for implants \varnothing 3.30 and 4.10SP have a M1.8 thread (like the thread of the implant platform). The transgingival caps for implants 3.75, 4.10 and 5.00 have a M2.0 thread.



implant diameter

**ø 3.30 mm and
ø 4.10SP mm**

**ø 3.75 mm and
ø 4.10 mm**

ø 5.00 mm


Transgingival healing caps
Anatomical emergence profile
Transgingival H. 2 mm



E-TMG-330-372



E-TMG-410-502




E-TMG-500-602


Transgingival healing caps
Anatomical emergence profile
Transgingival H. 3 mm



E-TMG-330-373




E-TMG-410-503




E-TMG-500-603


Transgingival healing caps
Anatomical emergence profile
Transgingival H. 5 mm



E-TMG-330-375



E-TMG-410-505



E-TMG-500-605

Transgingival healing caps
Straight emergence profile
Transgingival H. 2 mm



E-TMG-330-2

-

-

Transgingival healing caps
Straight emergence profile
Transgingival H. 3 mm



E-TMG-330-3

-

-

Transgingival healing caps
Straight emergence profile
Transgingival H. 5 mm



E-TMG-330-5

-

-

Impression and model phase

The components for the phase of taking the impression and creating the model are produced by the same maximum precision CNC machines that produce the respective implants; this ensures a real guarantee of precision from the point of view of tolerance and fidelity in the reproduction of the clinical situation. The open tray and closed tray transfers are made of Gr. 5 titanium, anodised according to the colour code of the reference implant platform.

Closed tray transfer: ideal for taking impressions with the closed tray technique, the repositioning face allows the transfer to be put back exactly in the impression material so as to reproduce the orientation of the connection precisely.

Open tray transfer: handy for use when there is a condition of parallelism between implants. The retentive design guarantees a steady grip in the upper portion of the impression.



..... **Analog:** anodised according to the colour code, like the transfers, to identify the different prosthetic platforms.



Open tray transfer




implant diameter

ø 3.30 mm and
ø 4.10SP mm


ø 3.75 mm and
ø 4.10 mm

ø 5.00 mm


Open tray transfer
Anatomical emergence
Connecting screw included



Open tray transfer
Straight emergence
Connecting screw included



Connecting screw
Supplied with the transfers, it can
also be ordered separately as a spare



WARNING: In case of prostheses on more than one post, it is recommended to glue the open tray transfers together with resin in order to guarantee stability and solidity of the impression.

Closed tray transfer




implant diameter

ø 3.30 mm and
ø 4.10SP mm


ø 3.75 mm and
ø 4.10 mm

ø 5.00 mm

Closed tray transfer
Connecting screw included



Connecting screw
Supplied with the transfers, it can
also be ordered separately as a spare



Analogues




implant diameter

ø 3.30 mm and
ø 4.10SP mm

ø 3.75 mm and
ø 4.10 mm

ø 5.00 mm

Analogue



Recommended tightening torque: 8-10 Ncm.

All measurements are given in mm, unless indicated otherwise

SIMPLE temporary posts

The SIMPLE prosthetic protocol calls for practical and simple solutions to create the temporary posts. The temporary posts can be used in a conventional way after the bone healing period, or immediately after surgical insertion of the implants, if conditions exist for immediate loading. Instead of using transgingival healing screws, depending on the prosthetic protocols adopted, it is possible to condition the tissues with the temporary prosthesis made using these posts.

PEEK is an extremely resistant and highly biodegradable polymer, easy to drill even chair-side. The titanium base, with an anatomical emergence profile, guarantees maximum connection precision. These posts have a connection hexagon which allows them to be repositioned; they are therefore ideal as a support for cement-retained single crowns.



In Simple aesthetic posts, the wider flaring of the transgingival profile, adaptable to any anatomy by drilling, simplifies the immediate aesthetic conditioning of the mucosa.

These posts do not have the non-rotational and repositioning connection hexagon; they are therefore useful for making multiple temporary structures to be screwed directly onto the implants.




implant diameter

**ø 3.30 mm and
ø 4.10SP mm**


**ø 3.75 mm and
ø 4.10 mm**

ø 5.00 mm


SIMPLE temporary posts in PEEK
With Gr. 5 titanium base
Anatomical emergence profile
Repositionable
With hexagon
Connecting screw included



E-MPSR-330




E-MPSR-410




E-MPSR-500


SIMPLE temporary posts in titanium
Anatomical emergence profile
Repositionable
With hexagon
Connecting screw included



E-MPSA-330-EX




E-MPSA-410-EX




E-MPSA-500-EX


SIMPLE temporary posts in titanium
Anatomical emergence profile
Non-repositionable
Without hexagon, rotating
Connecting screw included



E-MPSA-330




E-MPSA-410



E-MPSA-500


SIMPLE aesthetic temporary posts
in titanium
Wide emergence profile
Non-repositionable
Without hexagon, rotating
Connecting screw included



E-MPS-330



E-MPS-410



E-MPS-500

Connecting screw
Supplied with the temporary posts,
it can also be ordered separately as
a spare



VM-180
VM-180-10



VM-200
VM-200-10

Use
VM-200

Single pack
Pack of 10 pieces

WARNING: it is recommended always to use test screws for the laboratory phases and to keep the new screw supplied for the final fastening in the mouth.

Pre-made posts

Made of Gr. 5 titanium, these posts are subjected to a controlled passivation process that changes their surface colour: the result is a characteristic golden pale yellow colour. This colour is obtained through an anodising process and, therefore, there is no type of coating, which allows combining the advantages of a highly biocompatible surface with prosthetic reconstructions that are very aesthetically pleasing.

The golden colour makes them ideal for prosthetic rehabilitations of particular aesthetic importance.

Posts with passing screw angled at 15°.



The internal hexagon guarantees the non-rotational aspect and precise repositioning of the post.



| implant diameter | ø 3.30 mm and ø 4.10SP mm | ø 3.75 mm and ø 4.10 mm | ø 5.00 mm |
|---|---|---|---|
| Pre-made straight posts Repositionable Anatomical emergence Transgingival H. 1 mm Connecting screw included | ø 3.30... ø 3.75... E-MD-330-371 | ø 4.10... ø 5.00... E-MD-410-501 | ø 5.00... ø 6.00... E-MD-500-601 |
| Pre-made straight posts Repositionable Anatomical emergence Transgingival H. 2 mm Connecting screw included | ø 3.30... ø 3.75... E-MD-330-372 | ø 4.10... ø 5.00... E-MD-410-502 | ø 5.00... ø 6.00... E-MD-500-602 |
| Pre-made straight posts Repositionable Anatomical emergence Transgingival H. 4 mm Connecting screw included | ø 3.30... ø 3.75... E-MD-330-374 | ø 4.10... ø 5.00... E-MD-410-504 | ø 5.00... ø 6.00... E-MD-500-604 |
| Pre-made straight posts Repositionable Straight emergence Transgingival H. 1 mm Connecting screw included | ø 3.30... E-MD-330-1 | - | - |
| Pre-made straight posts Repositionable Straight emergence Transgingival H. 2 mm Connecting screw included | ø 3.30... E-MD-330-2 | - | - |
| Pre-made straight posts Repositionable Straight emergence Transgingival H. 4 mm Connecting screw included | ø 3.30... E-MD-330-4 | - | - |
| Pre-made posts angled at 15° Repositionable Anatomical emergence Connecting screw included | ø 3.30... ø 3.75... E-MAR15-330 | ø 4.10... ø 4.50... E-MAR15-410 | ø 5.00... ø 5.50... E-MAR15-500 |
| Pre-made posts angled at 15° Repositionable Straight emergence Connecting screw included | ø 3.30... E-MA15-330 | - | - |
| Connecting screw supplied with the pre-made posts, it can also be ordered separately as a spare Single pack Pack of 10 pieces | VM-180 VM-180-10 | VM-200 VM-200-10 | Use VM-200 |

Recommended tightening torque: 20-25 Ncm.

All measurements are given in mm, unless indicated otherwise.


Standard and Simple millable posts

These posts are made of grade 5 titanium and are able to respond to complex anatomical requirements in terms of both narrow prosthetic spaces and disparallel implants. Alongside the traditional prosthetic protocols, Sweden & Martina has developed innovative solutions in collaboration with distinguished professionals and universities. Among these, the SIMPLE Technique allows perfect conditioning of the mucosa starting with temporary posts (see page 50) and uses a millable post with large dimensions for making a primarily custom-built final prosthesis

SIMPLE millable posts have a very wide emergence profile which can be adapted to any anatomy obtained with SIMPLE temporary aesthetic posts in the immediate conditioning phase.



Straight millable posts, with their characteristic profile in the shape of an inverted cone, are suitable for angles of up to 10° and for small profiles.



Pre-unloaded posts allow very pronounced angles to be reached, up to 25°, limiting the time needed to reduce them by milling.






implant diameter

**ø 3.30 mm and
ø 4.10SP mm**


**ø 3.75 mm and
ø 4.10 mm**

ø 5.00 mm


Straight millable posts
Repositionable
With hexagon
Connecting screw included



E-MFD-330-50




E-MFD-410-60




E-MFD-500-75


Pre-unloaded millable posts
Repositionable
With hexagon
Connecting screw included



E-MFP-330-50



E-MFP-410-60



E-MFP-500-75

SIMPLE millable posts
Repositionable
With hexagon
Connecting screw included



E-MFS-330



E-MFS-410



E-MFS-500

Connecting screw
Supplied with the millable posts, it
can also be ordered separately as a
spare



**VM-180
VM-180-10**



**VM-200
VM-200-10**

Use
VM-200

Single pack
Pack of 10 pieces

WARNING: it is recommended always to use test screws for the laboratory phases and to keep the new screw supplied for the final fastening in the mouth.


Castable posts with alloy base

This solutions combines the simplicity of castable solutions with a gold alloy base that is highly biocompatible and highly resistant to corrosion. The melting point of the alloy is such as to preserve the base against dimensional alterations at the time of overcasting the castable part.



The non repositionable version presents an internal coupling cylinder which guarantees the simple insertion of multiple structures.

The repositionable version is suitable for creating single elements.




implant diameter
**ø 3.30 mm and
ø 4.10SP mm**
**ø 3.75 mm and
ø 4.10 mm**
ø 5.00 mm

Castable posts
With premade alloy base for
overcasting
Repositionable
With hexagon
Connecting screw included


E-UC-330-EX

E-UC-410-EX

E-UC-500-EX


Castable posts
With premade alloy base for
overcasting
Non-repositionable
Rotating
Connecting screw included


E-UC-330-ROT

E-UC-410-ROT

E-UC-500-ROT

Spare castable sleeves for castable
posts with alloy bases


E-UCC-330

E-UCC-410

E-UCC-500

Connecting screw
Supplied with the millable posts, it
can also be ordered separately as a
spare


**VM-180
VM-180-10**

**VM-200
VM-200-10**

Use
VM-200

Single pack
Pack of 10 pieces

WARNING: it is recommended always to use test screws for the laboratory phases and to keep the new screw supplied for the final fastening in the mouth.

For the technical specifications of the gold alloy "1" refer to page 89.

Recommended tightening torque: 20-25 Ncm.











All measurements are given in mm, unless indicated otherwise.

Fully castable posts

Sweden & Martina fully castable posts are made of PMMA, a resin that does not leave any residue during the casting process; like all the prosthetic components, these posts are made by turning too, with consequent respect of micrometric tolerances.





| implant diameter | ø 3.30 mm and ø 4.10SP mm | ø 3.75 mm and ø 4.10 mm | ø 5.00 mm |
|---|--|---|--|
| Fully castable posts Repositionable With hexagon Anatomical emergence Connecting screw included |  E-CCR-330-EX |  E-CCR-410-EX |  E-CCR-500-EX |
| Fully castable posts Repositionable With hexagon Straight emergence Connecting screw included |  E-CC-330-EX | - | - |
| Fully castable posts Non-repositionable Rotating Anatomical emergence Connecting screw included |  E-CCR-330-ROT |  E-CCR-410-ROT |  E-CCR-500-ROT |
| Fully castable posts Non-repositionable Rotating Straight emergence Connecting screw included |  E-CC-330-ROT | - | - |
| Connecting screw Supplied with the millable posts, it can also be ordered separately as a spare |  |  | Use VM-200 |
| Single pack Pack of 10 pieces | VM-180 VM-180-10 | VM-200 VM-200-10 | |

WARNING: it is recommended always to use test screws for the laboratory phases and to keep the new screw supplied for the final fastening in the mouth.

Prosthesis on intermediate abutments

These abutments are composed of a titanium base, which screws directly onto the implants, characterised by a small upper cone with a height of 0.70 mm, the same for all the connection diameters, which allows easy insertion and removal of the over-structures, even in case of slight disparallelism. Castable sleeves are supplied with the abutment, to be used for modelling and casting the over-structure, and the respective connecting screw. A small hexagon is found at the base of the cone which is used as a "key for screwing" the piece to the implant. To transfer the abutment into the oral cavity each package contains a practical plastic carrier (code AVV-ABUT-DG, not available separately). For the final fixing of the abutments to the implants, instead, use the special key, code AVV2-ABUT, not included in the surgical kits, but in the Screw-Kit of prosthetic drivers, and it can also be ordered separately (see table on page 37).

The torque recommended for fastening the direct screw-retained abutments is 25-30 Ncm.





| implant diameter | ø 3.30 mm and ø 4.10SP mm | ø 3.75 mm and ø 4.10 mm | ø 5.00 mm |
|---|------------------------------|----------------------------|-----------------------|
| Standard Abutments Non-repositionable with straight screwing Complete with castable sleeve and respective connecting screw Emergence H. 2.5 mm | | | Use E-ABUT-410-2.5 |
| Standard Abutments Non-repositionable with straight screwing Complete with castable sleeve and respective connecting screw Emergence H. 3.5 mm | | | Use E-ABUT-410-3.5 |
| Protection cap for abutments | | Use E-ABUT-VT | Use E-ABUT-VT |
| Spare castable sleeve, for abutments | | Use E-ABUT-CC | Use E-ABUT-CC |
| Universal spare prosthetic screws for fastening the prosthesis to the abutments | | Use VABUT | Use VABUT |
| Open tray transfer without hexagon, for abutments Connecting screw included | | Use E-TRABUT | Use E-TRABUT |
| Spare screws for open tray transfer of the abutments | | Use VTRABUT | Use VTRABUT |
| Abutment analogs | | Use E-ANABUT | Use E-ANABUT |

WARNING: it is recommended always to use test screws for the laboratory phases and to keep the new screw supplied for the final fastening in the mouth.

Recommended tightening torque for securing the prosthetic screws: 20-25 Ncm.

P.A.D. Disparallel Screwed Prosthesis

The P.A.D. systematics (Disparallel Screwed Prosthesis) was designed to facilitate the production of multiple screwed prostheses, even in the presence of very divergent implants and disparallel prosthetic emergences. The P.A.D. angled abutments, in particular, are the simplest and most reliable solution for implants positioned in distal saddles with high inclination. The P.A.D. prosthetic system is very versatile, starting from the wide range of straight abutments (available in various transgingival heights of 1.5, 3 and 4 mm), angled abutments (available with angles of 30° and 17° and transgingival heights of 3 and 5 mm), and a complete line of components necessary for producing the over-structures (transfers, analogues, sleeves, etc.).

P.A.D. posts: these all have the same upper cone, angled at 15°, which facilitates the manoeuvres of inserting and removing multiple screwed prostheses.



The upper cone allows further repositioning of the prosthetic structure by 15° on each side, which in the case of angled P.A.D. abutments are added to the angle of 17° or 30°. This characteristic allows easy management of disparallelism of up to 45° on each side.



| implant diameter | ø 3.30 mm and ø 4.10SP mm | ø 3.75 mm and ø 4.10 mm | ø 5.00 mm |
|---|-----------------------------|-----------------------------|------------------------|
| Straight P.A.D. abutments Direct screw-retained Transgingival H. 1.5 mm | E-PAD-AD330-15 | E-PAD-AD410-15 | Use E-PAD-AD410-15 |
| Straight P.A.D. abutments Direct screw-retained Transgingival H. 3 mm | E-PAD-AD330-30 | E-PAD-AD410-30 | Use E-PAD-AD410-30 |
| Straight P.A.D. abutments Direct screw-retained Transgingival H. 4 mm | - | E-PAD-AD410-40 | Use E-PAD-AD410-40 |
| P.A.D. abutment angled at 17° Transgingival H. 3 mm Connecting screw included | E-PAD-AA330-173 | E-PAD-AA410-173 | Use E-PAD-AA410-173 |
| P.A.D. abutment angled at 17° Transgingival H. 5 mm Connecting screw included | - | E-PAD-AA410-175 | Use E-PAD-AA410-175 |
| P.A.D. abutment angled at 30° Transgingival H. 3 mm Connecting screw included | E-PAD-AA330-303 | E-PAD-AA410-303 | Use E-PAD-AA410-303 |
| P.A.D. abutment angled at 30° Transgingival H. 5 mm Connecting screw included | - | E-PAD-AA410-305 | Use E-PAD-AA410-305 |
| Spare screws for fastening the P.A.D. angled abutments to the implants Supplied with the P.A.D., they can also be ordered separately | | | Use PAD-VM-200 |
| Single pack Pack of 10 pieces | PAD-VM-180 PAD-VM-180-10 | PAD-VM-200 PAD-VM-200-10 | |










To transfer **straight** abutments into the oral cavity each package contains a practical plastic carrier (code AVV-ABUT-DG, not available separately). For fixing the abutments to the implants, instead, use the special driver, code AVV2-ABUT. The torque recommended for fastening the direct screw-retained abutments is 25-30 Ncm.

To transfer **angled** abutments into the oral cavity, a special instrument is available (code PAD-CAR-ABUA), present in the Screw Kit and which can be ordered separately (see page 37). To stabilise the abutment in the instrument a screw must be used (code PAD-VTRAL-140), to be ordered separately. The recommended tightening torque is 20-25 Ncm.

WARNING: it is recommended always to use test screws for the laboratory phases and to keep the new screw supplied for the final fastening in the mouth.

All measurements are given in mm, unless indicated otherwise.





PROSTHETIC COMPONENTS

| code | description |
|--|---|
|  <p>PAD-CG</p> | <p>Protection cap for P.A.D. abutments in Gr. 5 titanium, to be used if the abutments are not immediately "provisionalised". Connecting screw included (code PAD-VP-140), available also as a spare, to be tightened at 8-10 Ncm</p> |
|  <p>PAD-CAP</p> | <p>Rotating caps in POM for direct impression taking on P.A.D. abutments</p> |
|  <p>PAD-CAP-EX</p> | <p>Non-rotating caps in POM for direct impression taking on P.A.D. abutments, with hexagon</p> |
|  <p>PAD-TRA</p> | <p>Open tray transfer in Gr. 5 titanium for P.A.D. abutments, rotating. Long transfer screw included (code PAD-VTRAL-140), suitable for open impression tray and available also as a spare</p> |
|  <p>PAD-TRA-EX</p> | <p>Open tray transfer in Gr. 5 titanium for P.A.D. abutments, with hexagon, non-rotating. Long transfer screw included (code PAD-VTRAL-140), suitable for open impression tray and available also as a spare</p> |
|  <p>PAD-VTRAL-140</p> | <p>Spare Spare screw for P.A.D. abutment transfer supplied with the transfers, it can be ordered separately as a spare</p> |
|  <p>PAD-ANA</p> | <p>Analog for P.A.D. abutment in Gr. 5 titanium. Long transfer screw included (code PAD-VTRAL-140), suitable for open impression tray and available also as a spare</p> |
|  <p>PAD-CC</p> | <p>Castable sleeves in PMMA for P.A.D. abutments, rotating. Connecting screw included (code PAD-VP-140), available also as a spare. Attention: The recommended tightening torque for all the over-structures obtained by casting to the abutments is 20-25 Ncm. However, before casting, care must be taken in the laboratory to ensure that the castable sleeves are not fastened to the models with a torque exceeding 8-10 Ncm, because polymers are not as resistant as metal</p> |
|  <p>PAD-CC-EX</p> | <p>Castable sleeves in PMMA for P.A.D. abutments, with hexagon, non-rotating. Connecting screw included (code PADVP- 140), available also as a spare. Attention: The recommended tightening torque for all the over-structures obtained by casting to the abutments is 20-25 Ncm. However, before casting, care must be taken in the laboratory to ensure that the castable sleeves are not fastened to the models with a torque exceeding 8-10 Ncm, because polymers are not as resistant as metal</p> |

Recommended tightening torque for securing the prosthetic screws: 20-25 Ncm.

All measurements are given in mm, unless indicated otherwise.






| code | description |
|--|---|
|  PAD-CP | PEEK sleeves for P.A.D. abutments, rotating. They are specifically for creating a temporary prosthesis or in cases when it is necessary to reline an old prosthesis for using as a temporary one. Connecting screw included (code PAD-VP-140), available also as a spare, to be tightened at 20-25 Ncm |
|  PAD-CP-EX | PEEK sleeves for P.A.D. abutments, with hexagon, non-rotating. They are specifically for creating a temporary prosthesis or in cases when it is necessary to reline an old prosthesis for using as a temporary one. Connecting screw included (code PADVP- 140), available also as a spare, to be tightened at 20-25 Ncm |
|  PAD-UC | Castable posts in PMMA with a pre-made base in gold alloy type "1", rotating, not repositionable, for overcasting on P.A.D. abutments. Connecting screw included (code PAD-VP-140), available also as a spare, to be tightened at 20-25 Ncm regardless of the working phase, because the head of the screw never rests on the PMMA, but always on the alloy base. The castable sleeve is also available as a spare (code A-CCUCR-330) |
|  PAD-VP-140 | Spare screw for P.A.D. abutment prosthetic components Also available in packs of 10 pieces (code PAD-VP-140-10) |

For the technical specifications of the gold alloy "1" and of PMMA, refer to pages 89 and 87 respectively.

Recommended tightening torque for securing the prosthetic screws: 20-25 Ncm.

P.A.D. components for relining and gluing technique

| code | description |
|--|---|
|  PAD-CT | Sleeves in Gr. 5 titanium for P.A.D. abutments, rotating. They are specifically for an immediate and final prosthetisation process or for relining an old prosthesis to be used as a temporary post. Connecting screw included (code PAD-VP-140), available also as a spare, to be tightened at 20-25 Ncm |
|  PAD-CT-EX | Sleeves in Gr. 5 titanium for P.A.D. abutments, with hexagon, non-rotating. They are specifically for an immediate and final prosthetisation process or for relining an old prosthesis to be used as a temporary post. Connecting screw included (code PAD-VP-140), available also as a spare, to be tightened at 25-30 Ncm |
|  PAD-CCEM | Castable posts in PMMA for cementing on titanium sleeves Effective for prosthetisation without residual tensions. |

P.A.D. prostheses for "D.P.F." Technique (Direct Prosthetic Framework)

P.A.D. abutments have proven to be a valid support for creating various simplified prosthetic protocols, including the creation of temporary posts for full-arch implant rehabilitations with immediate loading with a very simple and safe procedure. The D.P.F. components have been specially developed for creating a castable resin structure directly in the oral cavity that is absolutely passive, not restricted by connection geometries and with the additional advantage of being made without errors due to the taking of the impression and the development of the model. The intra-oral cementing of the metal truss obtained subsequently by casting allows the times for inserting the reinforced temporary post to be reduced to 8 hours after the end of surgery, while still maintaining the important properties of resistance and passivity during the first phase of implant loading. The temporary post created in this way can also be used as a positioning stent for making the final prosthesis.

STEP 1

The castable bar is fixed to the caps with a photo-polymerisable resin.



STEP 3








The product is cast in the laboratory, obtaining the metal structure for making the finished truss.



STEP 2

After photo-polymerisation of the truss, the product is removed from the oral cavity.



| code | description |
|---|---|
|  <p data-bbox="256 555 325 584">PAD-LV</p> | <p data-bbox="483 439 1465 555">Complete pack of all the prosthetic components for the "D.P.F." technique on a single P.A.D. abutment. The pack includes the titanium sleeve (PAD-CT-LV), the castable centring device (PAD-CC-LV), the anti-escape plug (PAD-TR-LV), the protective o-ring (PAD-ORING-LV) and the connecting screw (PAD-VP-140) to be tightened at 20-25 Ncm, available also as a spare.</p> |
|  <p data-bbox="244 741 339 770">PAD-CT-LV</p> | <p data-bbox="759 656 1190 712">Spare titanium sleeve for the "D.P.F." technique. The pack does not include the connecting screw.</p> |
|  <p data-bbox="244 931 339 958">PAD-CC-LV</p> | <p data-bbox="719 857 1230 887">Spare castable centring device for the "D.P.F." technique.</p> |
|  <p data-bbox="244 1153 339 1180">PAD-TR-LV</p> | <p data-bbox="751 1061 1198 1090">Spare anti-escape plug for the "D.P.F." technique.</p> |
|  <p data-bbox="225 1368 357 1397">PAD-ORING-LV</p> | <p data-bbox="799 1279 1150 1308">Spare o-ring for the "D.P.F." technique.</p> |
|  <p data-bbox="236 1585 344 1612">PAD-VP-140</p> | <p data-bbox="715 1485 1238 1541">Spare screw for P.A.D. abutment prosthetic components Also available in packs of 10 pieces (code PAD-VP-140-10)</p> |
|  <p data-bbox="264 1805 320 1832">BARC</p> | <p data-bbox="834 1715 1118 1744">Castable bar, L. 5 cm, Ø 2.2 mm</p> |

Individualised prosthesis Echo2

In case of an individualised prosthesis, maximum aesthetics and design flexibility are obtained with ECHO individual posts and with ECHO Direct Bridges, designed using the CAD-CAM technique and produced at the Sweden & Martina ECHO milling centre. The individual posts in titanium are an evolution of the standard millable posts and allow the prosthesis to adapt perfectly to the patient's gingival anatomy, which is difficult to obtain using traditional laboratory techniques. Posts completely made of zirconium oxide are also available, which are the most advanced individualised solution today since they present many advantages in terms of translucence of restoration, extreme customisation of the product, biocompatibility and absence of corrosion in the oral cavity, maximum connection precision, excellent resistance to occlusal loads, less invasiveness thanks to the perfect adaptation to the tissues and reduced session times.

All the individual components (posts, Direct Bridges, etc.) are sold complete with the necessary connecting screws. For zirconium components a washer is also supplied, made of a special high-resistance polymer to be placed between the head of the screw and its stop, with the function of absorbing and distributing forces between the parts. Screws and Peek washers are also available as spares.

The aluminium Scan-transfer is available for scanning (code E-CAMETRA330 and E-CAMETRA410, details on page 70).

Individual posts:

- titanium
- zirconium
- cobalt chrome



Screw-retained bar structures:


- milled cobalt chrome
- milled bio-titanium





Caps and bridges for prosthesis cemented on abutments:










- zirconium
- milled cobalt chrome
- sintered cobalt chrome
- biomedical resin for temporary posts
- PMMA
- fibreglass
- milled bio-titanium



Screw-retained bridge structures and DIRECT BRIDGES

- zirconium
- milled cobalt chrome
- milled bio-titanium
- PMMA
- fibreglass

PROSTHETIC COMPONENTS

| implant diameter | ø 3.30 mm and ø 4.10SP mm | ø 3.75 mm and ø 4.10 mm | ø 5.00 mm |
|--|---|--|--|
| ECHO Scan Transfer for implants (in aluminium ERGAL 7075) Connecting screw included |  E-CAMETRA330 |  E-CAMETRA410 |  Utilizzare E-CAMETRA410 |
| Spare screws for ECHO Scan Transfer for implants |  VM-180 VM-180-10 |  VM-200 VM-200-10 | Use VM-200 |
| Single pack Pack of 10 pieces | | | |
| Spare screws for fastening ECHO individual posts and prosthetic over- structures in zirconium oxide directly on implants (in Gr. 5 titanium, complete with lock ring washer) | - |  E-CAMTVABU200 | - |
| Spare screws for fastening ECHO individual posts and prosthetic over- structures in titanium and cobalt chrome directly on implants | Use VM-180 | Use VM-200 | Use VM-200 |
| Spare screws for fastening ECHO prosthetic over-structures in zirconium oxide directly on P.A.D. abutments (in Gr. 5 titanium, complete with lock ring washer) |  M1.4 PAD-VCAM-140 | Use PAD-VCAM-140 | Use PAD-VCAM-140 |
| Spare screws for fastening ECHO prosthetic over-structures in titanium and cobalt chrome on P.A.D. abutments |  M1.4 PAD-VP-140 | Use PAD-VP-140 | Use PAD-VP-140 |
| Spare lock ring washers for the head of the spare connecting screw, for ECHO individual posts in zirconium (in Classic Peek). |  | Use CAMPRON205-10 | Use CAMPRON205-10 |
| Pack of 10 pieces | CAMPRON205-10 | | |

WARNING: In case of prostheses on more than one post, it is recommended to glue the transfers together with resin in order to guarantee stability and solidity of the impression.

Recommended tightening torque for securing the prosthetic screws: 20-25 Ncm.



T-Connect

T-Connect supports can be used for making individual prostheses in zirconium with open CAD-CAM systems, including Echo2 by Sweden & Martina, without sacrificing micrometric precision in the join between the platforms that can be obtained with traditional components. Users of Echo2 can also choose to use T-Connect supports: the zirconium posts obtained in this way require a small support base in titanium which avoids contact between the zirconium oxide body and the implant platform.

For further information on the compatible systems, contact the CAD-CAM product specialist at Sweden & Martina.



T-Connect: allow a perfect join to be obtained between prosthesis and implant with zirconium structures milled and sintered in the laboratory.

implant diameter

ø 3.30 mm and
ø 4.10SP mm

ø 3.75 mm and
ø 4.10 mm

ø 5.00 mm

T-Connect in Gr. 5 titanium for individualised posts in zirconium
Connecting screw included



E-BASTZR-S-330



E-BASTZR-S-410



E-BASTZR-S-500

Connecting screw supplied with the T-Connect, it can also be ordered separately as a spare



VM-180
VM-180-10



VM-200
VM-200-10

Use VM-200

Single pack
Pack of 10 pieces

WARNING: it is recommended always to use test screws for the laboratory phases and to keep the new screw supplied for the final fastening in the mouth..

Recommended tightening torque for securing the prosthetic screws: 20-25 Ncm.

Locator abutment

Locator abutments are a patented and versatile prosthetic solution for attaching overdentures to dental implants easily and safely. The Locator system allows easily correcting misalignment of divergent implants by up to 40° (20° for each implant) in limited occlusal spaces; given the limited amount of space occupied, Locator abutments are perfect for patients with a removable prosthesis.








These abutments are made of Gr. 5 titanium and are available in different transgingival heights. The Locators are tightened at a torque of 25-30 Ncm using the special driver to be ordered separately (code 8926-SW, short, and code 8927-SW, long).






Retainers with different flexibilities are available for easily inserting into the steel cap with the use of a practical instrument for retainer insertion. The different retentive capacities are easily identified based on a colour code. The blue, pink and transparent retainers can be used on implants with angles of up to 10°; orange, red and green retainers are used on implants with an angle between 10° and 20°.

PROSTHETIC COMPONENTS





| implant diameter | ø 3.30 mm and ø 4.10SP mm | ø 3.75 mm and ø 4.10 mm | ø 5.00 mm |
|--|---|---|-----------|
| Locator abutment* Transgingival H. 1 mm |  1773 |  1741 | Use 1741 |
| Locator abutment* Transgingival H. 2 mm |  1774 |  1742 | Use 1742 |
| Locator abutment* Transgingival H. 3 mm |  1775 |  1743 | Use 1743 |
| Locator abutment* Transgingival H. 4 mm | - |  1744 | Use 1744 |

* *Locator abutments are medical devices manufactured and patented by Zest Anchors, Inc., 2061 Wineridge Place, Escondido, CA 92029, USA. Locator is a registered trademark of Zest Anchors, Inc. The European Agent for the purposes of MDD 93/42/EEC is Ventura Implant and Attachment Systems, 69 The Avenue, Ealing, London W13 8JR, England.*

Impression taking












| code | description |
|---|---|
|  8530 | Pack of 4 aluminium analogues for Locator abutments, one size for all platforms |



| code | description |
|---|---|
|  | Pack of 4 aluminium transfers for Locator abutments, one size for all platforms, and 4 black polyethylene retainers (LDPE 993I) with low retention included (code 8515). The retainers are also available as spares |
| 8505 | |
|  | Pack of 4 black polyethylene retainers (LDPE 993I) with low retention capacity for impression taking |
| 8515 | |
|  | Pack of 4 black nylon parallelism pins (LDPE 993I) for Locator abutments |
| 8517 | |
|  | Steel plate AISI 316L for measuring angles |
| 9530 | |





PROSTHETIC COMPONENTS

Plastic caps and retainers for Locator abutments*

| code | description |
|---|---|
|  <p>8519-2</p> | <p>Kit containing 2 Gr. 5 titanium caps, 2 spacer rings in silicon rubber, 2 black polyethylene retainers (LDPE 993I) with low retention capacity for impression taking and 2 nylon retainers for each of the 4 different retention capacities, designed for slight disparallelism (up to 10°)</p> |
|  <p>8540-2</p> | <p>Kit containing 2 Gr. 5 titanium caps, 2 spacer rings in silicon rubber, 2 black polyethylene retainers (LDPE 993I) with low retention capacity for impression taking and 2 nylon retainers for each of the 4 different retention capacities, designed for severe disparallelism (between 10° and 20°)</p> |
|  <p>8550-2</p> | <p>Kit containing 2 steel caps, 2 spacer rings in silicon rubber, 2 black polyethylene retainers (LDPE 993I) with low retention capacity for impression taking and 2 nylon retainers for each of the 4 different retention capacities, designed for slight disparallelism (up to 10°)</p> |
|  <p>8514</p> | <p>Pack of 20 spacer rings in silicon rubber, for the prosthesis relining phase</p> |
|  <p>8515</p> | <p>Pack of 4 black polyethylene retainers (LDPE 993I) with low retention capacity for impression taking</p> |
|  <p>8524</p> | <p>Pack of 4 transparent nylon retainers, retention 5 lb corresponding to 2268 g for slight disparallelism (up to 10°)</p> |
|  <p>8527</p> | <p>Pack of 4 pink nylon retainers, retention 3 lb corresponding to 1361 g for slight disparallelism (up to 10°)</p> |
|  <p>8529</p> | <p>Pack of 4 blue nylon retainers, retention 1.5 lb corresponding to 680 g for slight disparallelism (up to 10°)</p> |
|  <p>8547</p> | <p>Pack of 4 green nylon retainers, retention 4 lb corresponding to 1814 g for severe disparallelism between 10° and 20°)</p> |
|  <p>8548</p> | <p>Pack of 4 red nylon retainers, retention 1 lb corresponding to 450 g for severe disparallelism between 10° and 20°)</p> |
|  <p>8915</p> | <p>Pack of 4 orange nylon retainers, retention 2 lb corresponding to 907 g for severe disparallelism between 10° and 20°)</p> |

* *Locator Abutments are medical devices manufactured and patented by Zest Anchors, Inc., 2061 Wineridge Place, Escondido, CA 92029, USA. Locator is a registered trademark of Zest Anchors, Inc. The European Agent for the purposes of MDD 93/42/EEC is Ventura Implant and Attachment Systems, 69 The Avenue, Ealing, London W13 8JR, England.*



| code | description |
|--|---|
|  | <p>Locator Core Tool. Steel instrument composed of:</p> <ul style="list-style-type: none"> - handle - driver (8390) for screwing Locator abutments. - push-rod (8397) for inserting retainers into the caps - spare (8394) for Locator Core Tool retainer |
| 8393 | |
|  | <p>Steel push rod for inserting the retainers in the caps. Not necessary for those who already own or order the complete Locator Core Tool separately</p> |
| 8397 | |
|  | <p>Steel driver for abutment screwing/unscrewing. Not necessary for those who already own or order the Locator Core Tool separately</p> |
| 8390 | |
|  | <p>Retention jacket for the driver (8390) for transferring the Locator abutments into the oral cavity</p> |
| 8394 | |

Overdentures anchored with ball attachments

These attachments have a small hexagon at the base of the ball for engaging the driver (N.B.: the driver code BASCC-EX is not included in the surgical kit, it is in the Screw Kit and may be ordered separately). This driver is compatible with the system's dynamometric ratchet.





| implant diameter | ø 3.30 mm and ø 4.10SP mm | ø 3.75 mm and ø 4.10 mm | ø 5.00 mm |
|--|---------------------------|-------------------------|-------------------|
| Ball attachment Transgingival H. 1 mm | E-AS-330-1 | E-AS-410-1 | Use E-AS-410-1 |
| Ball attachment Transgingival H. 2 mm | E-AS-330-2 | E-AS-410-2 | Use E-AS-410-2 |
| Ball attachment Transgingival H. 4 mm | E-AS-330-4 | E-AS-410-4 | Use E-AS-410-4 |
| Analog of the ball attachment | ANAS | Use ANAS | Use ANAS |

code

description



BASCC-EX



Steel driver for ball attachments, with connector for dynamometric ratchet or digital connector.
(Not included in the surgical kit, included in the Screw Kit and can be ordered separately.)

Recommended tightening torque: 25-30 Ncm.






PROSTHETIC COMPONENTS

Accessories for overdentures on ball attachments


Polyamide caps for ball attachments

| code | description |
|---|--|
|  CAP-TFL-1 | Polyamide cap for ball attachments Ø 2.20 mm |
|  CONT-CAP-TFL-1 | Steel container for polyamide cap with outer Ø 4.80 mm. The total height is 3.20 mm |

Titanium caps for ball attachments





| code | description |
|---|---|
|  CAP-TIT-1 | Gr. 5 titanium cap complete with cap in two parts, titanium retention spring, and plastic mounting ring for ball attachments Ø 2.20 mm. The total height is 3.20 mm |
|  AN-CAP-TIT-1 | Spare plastic ring for titanium cap height 2.20 mm |
|  MOL1-CAP-TIT-1 | Spare retention spring for titanium caps, average hardness, steel, Ø 3.20 mm |
|  MOL2-CAP-TIT-1 | Spare retention spring for titanium cap, soft, for progressive adaptation of the prosthesis, steel, Ø 3.20 mm |
|  AVV-CAP-TIT-1 | Instrument for assembling and maintaining the titanium CAP-TIT-1 |

Caps in gold alloy for ball attachments

| code | description |
|--|---|
|  CAP-1 | Cap in gold alloy, complete with plastic positioning ring for ball attachments Ø 2.20 mm. The total height is 3.00 mm, and the outside diameter is 3.50 mm |



O-ring retention devices for ball attachments

| code | description |
|---|---|
|  99-440044* | Metal container in the shape of a ring for rubber O-rings. For ball attachments Ø 2.20 mm. The total height is 1.50 mm, and the outside diameter is 4.50 mm. Pack of 6 pieces |
|  99-443034* | Red ring in silicon for laboratory use, outside Ø 4.50 mm, height 1.50 mm. Pack of 12 pieces |
|  99-443035* | White ring in natural rubber, soft, outside Ø 4.50 mm, height 1.50 mm. Pack of 12 pieces |
|  99-443036* | Black ring in natural rubber, hard, outside Ø 4.50 mm, height 1.50 mm. Pack of 12 pieces |

* The retention O-ring for ball attachments are manufactured by Implant Direct Sybron International, 27030 Malibù Hills Road, Calabasas Hills, 91301 U.S.A. The European Agent for the purposes of MDD 93/42/EEC is Kerr Italia S.r.l., via Passanti 332, 84018 Scafati (SA) Italy.

Overdentures on bars

| code | description |
|---|---|
|  BARC-CAV-TIT | Castable bar, length 5.00 cm, height 3.00 mm, thickness 2.20 mm Ovoid-shaped profile with spacer |
|  CAV-TIT | Divisible bar attachment in titanium for oval bars height 3.00 mm, thickness 2.20 mm |
|  BARC | Castable bar, length 5.00 cm, Ø 2.20 mm |
|  CAV-375 | Bar attachment in gold alloy, for round bars with Ø 2.20 mm |

SFI-Bar


The total removable prosthesis supported by implants is acquiring more and more importance in prosthetic restorations. To enable patients to enjoy a substantially better quality of life on a long-term basis thanks to optimum retention, the principle of treatment with bars on implants has been radically reviewed, creating the SFI-Bar. The result is an exclusive solution that offers numerous clinical and economic advantages for patients, dentists and dental technicians in comparison with conventional bars.

The abutment driver for the SFI-Bar can be ordered separately, with code 07000046-SW. This driver is compatible with the system's dynamometric ratchet.

See pages 84-85 for the list of accessories available.

Chair-side and lab-side processing: Thanks to the extremely reduced production times, the SFI-Bar can be processed either chair-side or in the laboratory. In both cases the exclusive bar system with prefabricated elements is convincing for its almost universal possibilities of use, with costs reduced to a minimum in comparison with the traditional bar solution.

The bar prostheses made with the SFI-Bar can be supported by 3, 4, 5 or 6 implants, they take up less space in the mouth and present characteristics of biocompatibility. In case of repair it is possible to replace the individual elements without any problems.



Immediate loading: SFI-Bar offers the ideal basis for carrying out the immediate loading treatment, which exploits the complete absence of stress that can be obtained with this prosthetic protocol. This improves acceptance of the prosthesis by the patient and creates the conditions for optimum osseointegration.















Passive Fit: The stress-free position of a bar treatment is a fundamental requirement for the long-term clinical success of a prosthesis supported by implants. SFI-Bar offers the ideal basis for carrying out the immediate loading treatment, which exploits the complete absence of stress that can be obtained with this prosthetic protocol.




| implant diameter | ø 3.30 mm and ø 4.10SP mm | ø 3.75 mm and ø 4.10 mm | ø 5.00 mm |
|---|------------------------------|--|-----------|
| Abutment for SFI-Bar Transgingival H. 4 mm | - |  | - |
| | | 05001182 | |
| Abutment for SFI-Bar Transgingival H. 5 mm | - |  | - |
| | | 05001183 | |
| Abutment for SFI-Bar Transgingival H. 6 mm | - |  | - |
| | | 05001184 | |

Recommended tightening torque: 25-30 Ncm.






OVERDENTURES

| code | description | material | pack |
|---|---|---------------------------------------|--|
|  05000337 | Prosthetization kit on 2 implants, including: 2 large ball connectors, 2 connecting screws, 1 tubular bar | Gr. 5 titanium | |
|  05000338 | Prosthetization kit on 4 implants, including: 2 large ball connectors, 2 small ball connectors, 2 hemispherical shells, 4 connecting screws, 3 tubular bars | Gr. 5 titanium | |
|  05000668 | Additional prosthetization kit, including: 1 large ball connector, 1 connecting screw, 2 tubular bars | Gr. 5 titanium | |
|  05000382 | Tubular bar L. 20 mm | Gr. 5 titanium | |
|  05000383 | Large ball connector | Gr. 5 titanium | |
|  05000384 | Small ball connector | Gr. 5 titanium | |
|  05000385 | Hemispherical shell | Gr. 5 titanium | |
|  05000386 | Connecting screw | Gr. 5 titanium | |
|  05000344 | Asymmetrical female L. 50 mm | gold alloy "2" | 1 piece |
|  05000358 | Female complete with sheaths L. 47.5 mm | grade 5 titanium and plastic material | 1 bar with 12 segments, 6 yellow sheaths 6 red sheaths |
|  05000387 | Female L. 47.5 mm | Gr. 5 titanium | 1 bar with 12 segments |
|  05000388 | Yellow retentive sheath – low retention | plastic material | 6 pieces |
|  05000389 | Red retentive sheath – medium retention | plastic material | 6 pieces |
|  05000390 | Green retentive sheath – high retention | plastic material | 6 pieces |

Instruments

| code | description |
|---|---|
|  | Instruments kit |
| 07000108 | |
| 052082 | Brass wire for resilience To be inserted between the bar and the female during polymerisation to ensure a vertical translation of the prosthesis |
| 07000107 | Transfer axis L. 26 mm |
| 07000106 | Tube guide |
| 07000100 | Guide for cutting tubular bars |
| 07000114 | Flathead screwdriver for screwing the abutments onto the implants |
| 07000115 | Allen screwdriver for screwing the connecting screws onto the abutments |
| 070221 | Thomas key |
| 07000036 | Positioning pin for sheaths |
| 070198 | Set of activators for Elitor females |
| 070201 | Macro deactivator for Elitor females |
| 070347 | Forceps |
| 07000111 | Implant planner |
| 08000101 | Discs for cutting tubular bars |

Sweden & Martina auxiliary instruments

| code | description |
|---|---|
|  | Steel driver for abutments for SFI-Bar compatible with the dynamometric ratchet |
| 07000046-SW | |
|  | Mandrel for cutting discs, HP attachment |
| 305/050HP | |
|  | Abrasive corundum tip, HP attachment |
| 733MARRHP | |
|  | Rubber polishing head, HP attachment |
| 9503HP | |
|  | Rubber polishing head, CA attachment |
| 9503CA | |

MATERIAL COMPOSITION

GRADE 4 TITANIUM (cold worked)*

| | Element | Maximum allowed values (%) | Tolerance |
|-----------------------|----------|----------------------------|--|
| Chemical composition: | Nitrogen | 0.05 | +/- 0.02 |
| | Carbon | 0.08 | +/- 0.02 |
| | Hydrogen | 0.015 | +/- 0.002 |
| | Iron | 0.50 | +/- 0.01 (%<0.25) +/- 0.15 (%>0.25) |
| | Oxygen | 0.40 | +/- 0.02 (%<0.20) +/- 0.03 (%>0.20) |
| | Titanium | remainder | - |

Mechanical properties*

| | Minimum allowed values |
|------------------------|------------------------------|
| Tensile stress: | 680 MPa (N/mm ²) |
| Yield strength (0.2%): | 520 MPa (N/mm ²) |
| Elongation at yield: | 15 % |
| Section reduction: | 25 % |

* This technical information complies with the express specifications of the regulations in force for the use of grade 4 titanium in implantology:

- ASTM F67-06: Standard Specification for unalloyed titanium, for surgical implant applications.
- ISO 5832-2:1999: Implant for surgery – Metallic materials – Part 2: Unalloyed titanium.

Please note: the use of bars obtained from cold processing, for the production of Sweden & Martina Spa implants, allows the exploitation of the mechanical characteristics of tensile strength and yield strength about 15% higher than those that can be obtained with a hot process (respectively 550 MPa and 483 MPa).

GRADE 5 TITANIUM

| | Element | Maximum allowed values (%) | Tolerance |
|-----------------------|-----------|----------------------------|-----------|
| Chemical composition: | Nitrogen | 0.05 | +/- 0.02 |
| | Carbon | 0.08 | +/- 0.02 |
| | Hydrogen | 0.012 | +/- 0.002 |
| | Iron | 0.25 | +/- 0.10 |
| | Oxygen | 0.13 | +/- 0.02 |
| | Aluminium | 0.50÷6.50 | +/- 0.40 |
| | Vanadium | 3.50÷4.50 | +/- 0.15 |
| | Titanium | remainder | - |

Mechanical properties*

| | Minimum allowed values |
|--|------------------------------|
| Tensile stress (for bar diameters up to 44.45 mm): | 860 MPa (N/mm ²) |
| Yield strength (0.2%): | 795 MPa (N/mm ²) |
| Elongation at yield: | 10 % |
| Section reduction: | 25 % |

* This technical information complies with the express specifications of the regulations in force for the use of grade 5 titanium in implantology:

- ASTM F136-11: Standard Specification for wrought Titanium-6Aluminum-4Vanadium ELI (Extra low Interstitial) Alloy for surgical implant applications;
- ISO 5832-3:1996: Implants for surgery – Metallic materials – Part 3: Wrought titanium 6-aluminium 4-vanadium alloy.



PMMA

| | |
|---|----------------------------|
| Chemical designation: | Polymethylmethacrylate |
| Colour: | Transparent |
| Physical and mechanical properties | |
| Density (DIN 53479) | 1.18 g/cm ³ |
| Yield strength (DIN 53454) | 110 N/mm ² |
| Elongation at yield (DIN 53455) | 5.5 % |
| Modulus of elasticity (DIN 53457) | 115 N/mm ² |
| Tangential elastic modulus at 10 Hz (DIN 53445) | 3300 N/mm ² |
| BRINELL hardness ball falling h961/30 (DIN 53456) | 1700 N/mm ² |
| BRINELL hardness if the ball falls (DIN 53456) | 200 N/mm ² |
| Thermal properties | |
| Coefficient of linear extension for 0...50C (DIN VDE 0304/01) | 70·10 ⁻⁶ · 1/°C |
| Thermal conductivity (DIN 52612) | 0.19 W/m °C |
| Oven temperature | ≈ 160 °C |
| Regaining temperature | >80 °C |
| Maximum service temperature long term | 78 °C |
| VICAT temperature, proceeding B (DIN 53460) | 115 °C |
| ISO 75 flection resistance 1.80N/mm ² (DIN 53461) | 105 °C |
| Heat resistance according martens (DIN 53458) | 95 °C |
| Chemical properties | |
| Water absorption in weight increase after 1 day immersion (DIN 53495) | 0.3 % |

POM


| | |
|---|------------------------------|
| Chemical designation | Polyoxymethylene (copolymer) |
| Colour | Natural opaque |
| Physical and mechanical properties | |
| Density (DIN 53479): | 1.41 g/cm ³ |
| Tensile strength (DIN EN ISO 527-2) | 67 Mpa |
| Tensile strenght at yield (DIN EN ISO 527-2) | 67 Mpa |
| Elongation at yield (DIN EN ISO 527-2): | 9% |
| Elongation at break (DIN EN ISO 527-2): | 32 % |
| Modulus of elasticity (flexural test) (DIN EN ISO 178): | 2800 Mpa |
| Modulus of elasticity (tensile test) (DIN EN ISO 527-2): | 2800 MPa |
| Ball indentation hardness (ISO 2039-1): | 165 MPa |
| Impact strength (Charpy) (DIN EN ISO 179-1eU): | Not broken |
| Compression modulus (EN ISO 604): | 2300 Mpa |
| Thermal properties | |
| Melting temperature (DIN 53765): | 166 °C |
| Glass transition temperature (DIN 53765): | -60 °C |
| Service temperature short term: | 140 °C |
| Service temperature long term: | 100 °C |
| Specific heat (ISO 22007-4): | 1.4 J/(gK) |
| Thermal conductivity (ISO 22007-4): | 0.39 W/ (mK) |
| Thermal expansion (CLTE) 23-60°C (DIN EN ISO 11359-1;2): | 13·10 ⁻⁵ /K |
| Thermal expansion (CLTE) 23-100°C (DIN EN ISO 11359-1;2): | 14·10 ⁻⁵ /K |
| Chemical properties | |
| Water absorption 24h / 96h (23°C) (DIN EN ISO 62) | 0.05/0.1% |

COMPOSIZIONE DEI MATERIALI

| PEEK *(tested on the same quantity of material) | Radiopaque | Classic |
|--|------------------------|-----------------------|
| Chemical designation | Polyetheretherketone | Polyetheretherketone |
| Colour | Cream white opaque | Cream white opaque |
| Physical and mechanical properties | | |
| Density: | 1.65 g/cm ³ | 1,4 g/cm ³ |
| Modulus of elasticity (tensile test) (DIN EN ISO 527-2): | 5200 MPa | 4100 MPa |
| Tensile strength (DIN EN ISO 527-2): | 77 MPa | 97 MPa |
| Tensile strength at yield (DIN EN ISO 527-2): | 77 MPa | 97 MPa |
| Elongation at yield (DIN EN ISO 527-2): | 2% | 5% |
| Elongation at break (DIN EN ISO 527-2): | 2% | 13% |
| Flexural strength (DIN EN ISO 178): | 178 MPa | 174 MPa |
| Modulus of elasticity (flexural test) (DIN EN ISO 178): | 5000 MPa | 4000 MPa |
| Compression modulus (EN ISO 604): | 4000 MPa | 3500 MPa |
| Thermal properties | | |
| Glass transition temperature (DIN 53765): | - | 150 °C |
| Service temperature short term: | 300 °C | 300 °C |
| Service temperature long term: | 260 °C | 260 °C |
| Chemical properties | | |
| Water absorption 24h / 96h (23°C) (DIN EN ISO 62): | - | 0.02/0.03 % |



| GOLD ALLOY | Gold alloy 1 | Gold alloy 2 | Gold alloy 3 |
|--|--|---|---|
| Chemical designation | Gold alloy 1 | Gold alloy 2 | Gold alloy 3 |
| Colour | White | Yellow | Yellow |
| Composition | | | |
| Au | 60 % | > 68.60 % | 70 % |
| Pt | 24 % | 2.45 % | 8.5 % |
| Pd | 15 % | 3.95 % | - |
| Ir | 1 % | 0.05 % | 0.10 % |
| Ag | - | 11.85 % | 13.40 % |
| Cu | - | 10.60 % | 7.50 % |
| Zn | - | 2.50 % | 0.50 % |
| Au + Pt Group Metals | - | 75.35 % | - |
| Ru | - | - | - |
| Physical and mechanical properties | | | |
| Density: | 18.1 g/cm ³ | 15.0 g/cm ³ | 15.7 g/cm ³ |
| Melting range: | 1400 ÷ 1460 °C | 880 ÷ 940 °C | 895 ÷ 1010 °C |
| Modulus of elasticity (tensile test): | 115 GPa | 97 GPa | 100 GPa |
| Vickers Hardness HV1 (Gold alloy 1) HV5 (Gold alloy 2, Gold alloy 3) | 160 (annealed) 250 (hardened) 220 (after deformation) 240 (after casting) | > 240 | 170 (annealed) 295 (hardened) 280 (selfhardening) |
| Proof stress Rp0.2 | 400 MPa (annealed) 700 (after deformation) 800 (after casting) | > 710 MPa (cold worked) 410 Mpa (soft) 680 Mpa (hardened) | 380 MPa (annealed) 730 (after deformation) |
| Yield strenght Rm | 600 Mpa (annealed) 850 (hardened) 850 (after deformation) | >790 MPa (cold worked) 535 MPa (soft) 780 MPa (hardened) | - |
| Elongation in % | 20 % (annealed) 15 (hardened) 1 (after deformation) | > 4 % (cold worked) 35% (soft) 12% (hardened) | - |



Rev. 11-14



www.sweden-martina.com

Sweden & Martina S.p.A.
Via Veneto, 10 - 35020 Due Carrare (PD), Italy
Tel. +39.049.9124300
Fax +39.049.9124290

info@sweden-martina.com

Sweden & Martina Mediterranea S.L.
Sorolla Center, Oficina 801
Av.da Cortes Valencianas 58, 8pl -46015-Valencia, España
Tel. +34.96.3525895

info.es@sweden-martina.com
Numero gratuito 900993963

Except when indicated in the internal pages, the implants, the prosthetic components and the surgical subject of this catalogue are Medical Devices manufactured by Sweden & Martina S.p.A. , they conform to the UNI EN ISO 9001:2008 / UNI EN 13485:2012 standards and are certified with the CE Mark (Class I) and CE 0476 mark (Class IIA and class IIB) in compliance with European Medical Device Directive No. 93/42 and European Directive No. 2007/47.

All contents are updated at the time of printing. Please check in the headquarters whether any updates are available.

The contents are up to date at the time of publication.
Check with the company for any subsequent updates.