# Silk

**Sutura chirurgica in seta naturale plurifilamento intrecciata e rivestita in silicone. Non assorbibile.**  
**Sutura quirúrgica de seda natural polifilamento trenzada y revestida en silicona. No absorbible.**

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### Descrizione

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## INSTRUCTION FOR USE

### Description

**SILK** suture is a sterile medical device. **SILK** is a natural suture of animal origin composed of fibroin (a protein deriving from species bombyx mori silkworm larvae). The suture **SILK** dyed black (Logwood Black C.I.75290) is obtained from the braiding of elementary multifilaments and is coated of silicone. **SILK** sutures are available in a range of gauge sizes (USP/EP) and lengths, attached to stainless steel needles of different sizes and shapes. The **SILK** suture is classified by the European Pharmacopoeia as non-absorbable, non-absorbable sutures permanently maintain the ability to mechanically support the edges of the wound for very long times. However, they can undergo a slow and gradual metabolization which in any case does not compromise their mechanical characteristics over time. For the **SILK** suture, a loss of tensile strength has been noted after about one year as a consequence of degradation «in vivo». As with any foreign body, it may cause slight inflammatory tissue re-action, followed by gradual encapsulation of the suture from fibrous connective tissue.

### Indications

**SILK** suture is indicated for use in any procedure where the use of a non-absorbable suture is required for tissue approximation and/or ligation.

### Contraindications

Due to the gradual loss tensile strength which may occur after prolonged periods "in vivo", **SILK** suture should not be used where permanent retention of the suture strength is required.

### Precautions

**SILK** sutures are to be used by qualified personnel experienced in technical and surgical procedures where the use of a non-absorbable suture is required. Care should be taken when handling the suture material to avoid any bending or flattening to the thread due to the application of surgical instruments such as forceps or needleholders. Greater knot security requires standard surgical technique of flat and square knots, as may be indicated by the circumstances and the surgeon's experience. Surgical sutures are disposable medical devices and must be disposed of as hospital waste materials.

### Adverse reaction

Adverse reactions associated with the use of this suture device include slight local inflammatory reactions as typical reaction to a foreign body, wound dehiscence, lack of adequate support to approximation of tissues subject to growth, stretching or extension unless additional support is used, lack of adequate support to tissue approximation in elderly, malnourished, weakened, anemic, obese or diabetic patients and in all conditions which may delay wound healing process.

### Content

36 sterile envelopes

### Warnings

Keep in cool dry place, away from any direct heat.  
 Expiration date is shown on each single envelope and on the outer box.  
 Disposable  
 Do not use if package is damaged.

### Sterility

**SILK** sutures are sterilized by Gamma Ray.  
 Sterile unless package is opened or damaged.  
 Do not resterilize.

Reserved for the exclusive use by surgical medical personnel, in surgery or in hospital setting.  
 Keep out of reach of children.

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**Ergon Sutramed Srl**
**Sweden & Martina S.p.A.**