



GUIDED SURGERY

ECHOPLAN PRO KIT


Fast and simple guided surgery




ECHOPLAN PRO KIT A: For guided insertion of Premium One and Prama cylindrical implants



ECHOPLAN PRO KIT C: For guided insertion of CSR implants



ECHOPLAN PRO KIT E: For guided insertion of Outlink² cylindrical implants with external hexagon



ECHOPLAN PRO KIT S: For guided insertion of tapered Shelta, Syra, Prama RF and CSR RF SL implants



Autoclavable tray for organising, sterilising and storing guided implant mounts



PRODUCT DESCRIPTION

The Sweden & Martina Echoplan PRO line includes four surgical kits with which it is possible to insert implants with cylindrical and conical morphologies that have either a Collex (internal hexagon and prosthetic support collar), a conical or an external hexagon connection. The tray's instrument path is intuitive and direct, for quick handling by the clinician or assistant.

ECHOPLAN PRO A: this kit includes all the instruments needed for guided preparation of implantation sites for Premium One and Prama implants, in hts. from 8.50 mm to 15.00 mm. There are also slots in the tray for burrs used for the preparation of short implants ht. 6.00 mm and ht. 7.00 mm as well as tall implants ht. 18 mm, which are an available option in supplementary sets or sold individually.

ECHOPLAN PRO C: this kit includes all the instruments needed for guided preparation of the implantation sites for CSR implants in hts. from 8.50 mm to 15.00 mm. There are also slots in the tray for burrs used for the preparation of short implants ht. 6.50 mm and of tall implants ht. 18.00 mm, which are an available option in supplementary sets or sold individually.

ECHOPLAN PRO E: this kit includes all the instruments needed for guided preparation of the implantation sites for Outlink² implants in hts. from 8.50 mm to 15.00 mm. The tray also has slots for burrs used for the preparation of short implants ht. 5.00 mm and 7.00 mm and of tall implants ht. 18.00 mm, which are an available option in supplementary sets or sold individually.

ECHOPLAN PRO S: this kit includes all the instruments needed for guided preparation of the implantation sites for Shelta, Syra, Prama RF and CSR RF SL implants in hts. from 8.50 mm to 15.00 mm. There are also slots in the tray for burrs used for the preparation of short implants ht. 4.00 mm, 5.00 mm, and 6.00 mm and of tall implants ht. 18.00 mm, which are an available option in supplementary sets or sold individually.
