


















Explanation of symbols on Sweden & Martina Labelling

		implant packaging	surgical instrument packs	prosthesis packaging
	Caution! See instruction for use	✓	✓	✓
	Batch number	✓	✓	✓
	Code	✓	✓	✓
	Manufacturer	✓	✓	✓
	Country of manufacture	✓	✓	✓
	UDI code, Unique Device Identification	✓	✓	✓
	Medical Device	✓	✓	✓
	Consult instruction for use www.sweden-martina.com	✓	✓	✓
	CE marking <i>Where applicable: The identification number of the Notified Body shall follow this symbol.</i>	✓	✓	✓
	American federal law restricts this device to sale by or by order of a professional practitioner	✓	✓	✓
	Do not re-sterilize	✓		
	Disposable product, do not reuse	✓		✓
	Do not use if the packaging is damaged	✓	✓	✓
	Non-sterile product		✓	✓
	Sterilized with ionizing radiation	✓		
	Single sterile barrier system with protective packaging inside	✓		
	Expiry date after which the product must not be used	✓		