

PREMIUM ONE SHELTA

Product overview





sweden & martina

Premium One implants



h	ø 3.30 ●	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
7.00	-	-	AS-ZT-425-070	AS-ZT-500-070
8.50	A-ZT-330-085	A-ZT-380-085	AS-ZT-425-085	AS-ZT-500-085
10.00	A-ZT-330-100	A-ZT-380-100	AS-ZT-425-100	AS-ZT-500-100
11.50	A-ZT-330-115	A-ZT-380-115	AS-ZT-425-115	AS-ZT-500-115
13.00	A-ZT-330-130	A-ZT-380-130	AS-ZT-425-130	AS-ZT-500-130
15.00	A-ZT-330-150	A-ZT-380-150	AS-ZT-425-150	AS-ZT-500-150
18.00	-	A-ZT-380-180	AS-ZT-425-180	-
Surgical cover screw				
	A-VT-330	A-VT-380	SH-VT-425-BL	SH-VT-500-VI

Shelta implants



h	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
8.50	SH-ZT-380-085	SH-ZT-425-085	SH-ZT-500-085
10.00	SH-ZT-380-100	SH-ZT-425-100	SH-ZT-500-100
11.50	SH-ZT-380-115	SH-ZT-425-115	SH-ZT-500-115
13.00	SH-ZT-380-130	SH-ZT-425-130	SH-ZT-500-130
15.00	SH-ZT-380-150	SH-ZT-425-150	SH-ZT-500-150
Surgical cover screw			
	SH-VT-380-VE	SH-VT-425-BL	SH-VT-500-VI

Shelta SL implants



h	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
8.50	SH-ZT-380SL-085	SH-ZT-425SL-085	SH-ZT-500SL-085
10.00	SH-ZT-380SL-100	SH-ZT-425SL-100	SH-ZT-500SL-100
11.50	SH-ZT-380SL-115	SH-ZT-425SL-115	SH-ZT-500SL-115
13.00	SH-ZT-380SL-130	SH-ZT-425SL-130	SH-ZT-500SL-130
15.00	SH-ZT-380SL-150	SH-ZT-425SL-150	SH-ZT-500SL-150
Surgical cover screw			
	SH-VT-380-VE	SH-VT-425-BL	SH-VT-500-VI

Surgical kits

Premium One surgical kit



ZPREMIUM-ONE-INT grommetless surgical kit including all the instruments for Premium One implants


For further information on the surgical kit and the surgical procedures, please consult the surgical manual MC-IMP-PREMIUM-ONE-E.

Drilling kit for short Premium One implants



ZSHORTY-INT Shorty drilling kit including for short Premium One implants


Shelta surgical kit



ZSHELTA-INT grommetless surgical kit including all the instruments for Shelta implants


For further information on the surgical kit and the surgical procedures, please consult the surgical manual MC-IMP-SHELTA-USA.

Premium One and Shelta surgical kit



ZPRESH-INT grommetless surgical kit including all the instruments for Premium One and Shelta implants


Premium One guided surgery surgical kit



ZGS-PRO-A-INT grommetless surgical kit including all the instruments for guided surgery insertion of Premium One implants

For further information on the surgical kit and the surgical procedures, please consult the surgical manual MC-UD-EPLANPROA-USA.

Shelta and Shelta SL guided surgery surgical kit



ZGS-PRO-S-INT grommetless surgical kit including all the instruments for guided surgery insertion of Shelta and Shelta SL implants

For further information on the surgical kit and the surgical procedures, please consult the surgical manual MC-UD-EPLANPROS-E.

Instruments included in Premium One surgical kit

Precision and intermediate drills

	FS-230	precision drill
	FG-200/280XS	intermediate drill 200/280
	FG-330/425XS	intermediate drill 330/425

Pilot drill and related stops

	FPT3-200-LXS	
	h	
	7.00	STOP4-200-070
	8.50	STOP4-200-085
	10.00	STOP4-200-100
	11.50	STOP4-200-115
	13.00	STOP4-200-130
	15.00	STOP4-200-150

Bone taps with hexagonal connector

	ø 3.30	ø 3.80	ø 4.25	ø 5.00
A-MS-330	A-MS-380	A-MS-425	A-MS-500	

Screwdrivers

	HSMXS-20-DG	digital, extra short
	HSM-20-DG	digital, short
	HSML-20-DG	digital, long
	HSM-20-EX	with hexagonal connector
	HSML-20-EX	with hexagonal connector, long
	HSM-20-CA	with right angle shank

Dynamometric ratchet

	CRIS-KIT
---	----------

X-ray templates



AS-L100	real dimensions
AS-L120	dimensions increased by 20%
AS-L130	dimensions increased by 30%

Premium One final drills and related stops

		ø 3.30 	ø 3.80 	ø 4.25 	ø 5.00 
		FFT3-280-LXS	FFT3-300-LXS	FFT3-340-LXS	FFT3-425-LXS
	h	ø 3.30 	ø 3.80 	ø 4.25 	ø 5.00 
	7.00	STOP4-280-070	STOP4-300-070	STOP4-340-070	STOP4-425-070
	8.50	STOP4-280-085	STOP4-300-085	STOP4-340-085	STOP4-425-085
	10.00	STOP4-280-100	STOP4-300-100	STOP4-340-100	STOP4-425-100
	11.50	STOP4-280-115	STOP4-300-115	STOP4-340-115	STOP4-425-115
	13.00	STOP4-280-130	STOP4-300-130	STOP4-340-130	STOP4-425-130
	15.00	STOP4-280-150	STOP4-300-150	STOP4-340-150	STOP4-425-150

Easy Insert Drivers for implant insertion

	EASYC4-EX230-CA	short
	EASYL4-EX230-CA	long
	EASY4-EX230-EX	with hexagonal connector

Extensions and adapters

	BPM-15	extension for instruments with hexagonal connector
	PROF-CAL3	extension for surgical drills
	AVV3-MAN-DG	digital knob
	B-AVV-CA3	mechanical adapter for instruments with hexagonal connector

Parallelism pin with depth laser marking

	PPTL-2-28
---	-----------

Instruments included in Shelta surgical kit

Precision drill



FS-230

precision drill

Pilot drill, intermediate drill and related stops



CSR-FP-200

pilot drill



SE-FK250

intermediate tapered drill

h

8.50 CSR-STOP-2028-085

10.00 CSR-STOP-2028-100

11.50 CSR-STOP-2028-115

13.00 CSR-STOP-2028-130

15.00 CSR-STOP-2028-150

Bone taps with right angle shank



ø 3.80

ø 4.25

ø 5.00

SH-MS-380-CA

SH-MS-425-CA

SH-MS-500-CA

Screwdrivers



HSMXS-20-DG

digital, extra short

HSM-20-DG

digital, short

HSML-20-DG

digital, long



HSM-20-EX

with hexagonal connector

HSML-20-EX

with hexagonal connector, long



HSM-20-CA

with right angle shank

Dynamometric ratchet



CRI5-KIT

Easy Insert Drivers for implant insertion



EASYC4-EX230-CA

short



EASYL4-EX230-CA

long



EASY4-EX230-EX

with hexagonal connector

Shelta final drills and related stops



h

ø 3.80

ø 4.25

ø 5.00

8.50 SH-FK380-085 SH-FK425-085 SH-FK500-085

10.00 SH-FK380-100 SH-FK425-100 SH-FK500-100

11.50 SH-FK380-115 SH-FK425-115 SH-FK500-115

13.00 SH-FK380-130 SH-FK425-130 SH-FK500-130

15.00 SH-FK380-150 SH-FK425-150 SH-FK500-150



ø 3.80

ø 4.25

ø 5.00

SH-STOP4-FK380 SH-STOP4-FK425 SH-STOP4-FK500

Replies



h

ø 3.80

ø 4.25

ø 5.00

8.50 SH-380-085-RP SH-425-085-RP SH-500-085-RP

10.00 SH-380-100-RP SH-425-100-RP SH-500-100-RP

11.50 SH-380-115-RP SH-425-115-RP SH-500-115-RP

13.00 SH-380-130-RP SH-425-130-RP SH-500-130-RP

15.00 SH-380-150-RP SH-425-150-RP SH-500-150-RP

Extensions and adapters



BPM-15

extension for instruments with hexagonal connector



PROF-CAL3

extension for surgical drills



AVV-CA-DG-EX


adapter for right angle shank instruments

Parallelism pins with depth laser marking



PPTL-2-28

X-ray templates



SH-L100

real dimensions

SH-L120

dimensions increased by 20%

SH-L130

dimensions increased by 30%

Optional instruments (not included in Premium One and Shelta surgical kit)


Premium One implants initial drill and related stops

	FFT3-250-LXS	cylindrical drill
h		
	7.00 STOP4-250-070	
	8.50 STOP4-250-085	
	10.00 STOP4-250-100	
	11.50 STOP4-250-115	
	13.00 STOP4-250-130	
	15.00 STOP4-250-150	

Intermediate drill

	FG-200/280XS**	ø 2.00 - ø 2.80
---	----------------	-----------------

Premium One short implants initial drill and related stops

	FPS-200	pilot drill
	STOP4-200-070	stop for pilot drill
	FGS-200/300	intermediate drill 200/300
	FGS-340/425	intermediate drill 340/425
	FGS-425/540	intermediate drill 425/500

Premium One short implants final drills and related stops

	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
	FFS-300	FFS-340	FFS-425
	STOP4-300-070	STOP4-340-070	STOP4-425-070

Countersink drills

	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
	A-FCC-380	A-FCC-425	A-FCC-500

Premium One implants osteotomes with flat tip

	E-OS-020-PP	ø tip 0.20
	E-OS-090-PP	ø tip 0.90
	E-OS-160-PP	ø tip 1.60
	E-OS-200-PP	ø tip 2.00
	E-OS-240-PP	ø tip 2.40


Shelta implants osteotomes with flat tip

	h	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
	8.50-10.00	SH-OS-380-100-PP	SH-OS-425-100-PP	SH-OS-500-100-PP
	11.50	SH-OS-380-115-PP	SH-OS-425-115-PP	SH-OS-500-115-PP
	13.00	SH-OS-380-130-PP	SH-OS-425-130-PP	SH-OS-500-130-PP
	15.00	SH-OS-380-150-PP	SH-OS-425-150-PP	SH-OS-500-150-PP

Shelta implants osteotomes with round tip

	h	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
	8.50-10.00	SH-OS-380-100-PR	SH-OS-425-100-PR	SH-OS-500-100-PR
	11.50	SH-OS-380-115-PR	SH-OS-425-115-PR	SH-OS-500-115-PR
	13.00	SH-OS-380-130-PR	SH-OS-425-130-PR	SH-OS-500-130-PR
	15.00	SH-OS-380-150-PR	SH-OS-425-150-PR	SH-OS-500-150-PR

Bone profilers

transgingival h 3.00				
	ø 3.30 ●	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
	A-PAD-PS330-L	A-PAD-PS380-L	A-PAD-PS425-L	A-PAD-PS500-L
transgingival h 5.00				
	ø 3.30 ●	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
	A-PAD-PS330-S	A-PAD-PS380-S	A-PAD-PS425-S	A-PAD-PS500-S
	A-PAD-GUI-PS-230	guide cylinder		

Drills for distal sectors

	ø	FPT5-200-LXS	ø	FFT5-330-LXS
	2.00		3.30	
	2.80	FFT5-280-LXS	3.40	FFT5-340-LXS
	2.90	FFT5-290-LXS	3.60	FFT5-360-LXS
	3.00	FFT5-300-LXS	4.25	FFT5-425-LXS
	3.20	FFT5-320-LXS	4.45	FFT5-445-LXS

Premium One implants bone taps

standard, with right angle shank				
	ø 3.30 ●	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
	A-MS-330-CA	A-MS-380-CA	A-MS-425-CA	A-MS-500-CA
short, with hexagonal connector				
	ø 3.30 ●	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
	A-MS-330	A-MS-380	A-MS-425	A-MS-500

Shelta implants bone taps with hexagonal connector

	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
	SH-MS-380	SH-MS-425	SH-MS-500

Extensions and adapters

	AVV-CA-DG-EX*	adapter for right angle shank instruments
	B-AVV-CA3**	mechanical adapter for instruments with hexagonal connector
	AVV3-MAN-DG**	digital knob

Mounter

	MOU-EX230	mounter
	CM2	mounter stop key

P.A.D. carrier

	PAD-CAR
--	---------

Depth gauge

	PROF3
---	-------

Torque wrench with control lever

	TWL
--	-----

Spare O-ring

	ORING180-088
---	--------------

Parallelism pins

	PP-2/28	without depth laser marking
	PPTS-2-28	with depth laser marking, small
	PPS-2/3	for Premium One short implants

Extraction tools for intraoral removing of implants

	BC-EX230	short
	BL-EX230	long

Other prosthetic screwdrivers

	HSMXL-20-EX	with hexagonal connector, extra long
	HSMXS-20-CA	with right angle shank, extra short
	AVV-CAP-TIT-1	instrument for maintaining of the titanium cap for ball attachments
	BASCC-EX	for ball attachments, with hexagonal connector
	AVV2-ABUT	for straight P.A.D. abutments
	8926-SW	for Locator abutment, short
	8927-SW	for Locator abutment, long

*Not included in Premium One surgical kit, included in Shelta surgical kit.
**Not included in Shelta surgical kit, included in Premium One surgical kit

Implant platforms and prosthetic combinations

The connection platforms of the various implant sizes of Premium One and Shelta are characterized by the presence of the same hexagon and external Collex One collar, which allows the choice of the size of the abutment. According to the chosen prosthetic protocol, a Platform Switching can be made or a matching size abutment can be used.

ø implant	3.30 	3.80 	4.25 	5.00 
ø compatible prosthetic component	3.30 	3.30 	3.30 	3.30 
	-	3.80 	3.80 	3.80 
	-	-	4.25 	4.25 
	-	-	-	5.00 

Please note

Prosthetic components with ø 3.30 mm determine Platform Switching with ø 3.80, 4.25 and 5.00 mm. These components should only be used for single crowns in frontal sectors (premolars excluded), and in distal sectors only to support multiple restorations.

Tissues conditioning

Healing abutments

Straight healing abutments			
	transgingival height	ø 3.30 	ø 3.80 
	2.00	A-TMG-330-2	A-TMG-380-2
	3.00	A-TMG-330-3	A-TMG-380-3
	5.00	A-TMG-330-5	A-TMG-380-5
	7.00	A-TMG-330-7	-

Healing abutments with anatomical emergence					
	transgingival height	ø 3.30 	ø 3.80 	ø 4.25 	ø 5.00 
	2.00	A-TMGR-330-2	A-TMGR-380-2	AS-TMGR-425-2	AS-TMGR-500-2
	3.00	A-TMGR-330-3	A-TMGR-380-3	AS-TMGR-425-3	AS-TMGR-500-3
	5.00	A-TMGR-330-5	A-TMGR-380-5	AS-TMGR-425-5	AS-TMGR-500-5
	7.00	-	A-TMGR-380-7	AS-TMGR-425-7	AS-TMGR-500-7

Healing abutments with coronal enlargement			
	transgingival height	ø 3.30 	ø 3.80 
	2.00	A-TMTCR-330-2	A-TMTCR-380-2
	3.00	A-TMTCR-330-3	A-TMTCR-380-3
	5.00	A-TMTCR-330-5	A-TMTCR-380-5

Healing abutments with wide profile			
	transgingival height	ø 3.30 	ø 3.80 
	2.00	A-TMGB-330-2	A-TMGB-380-2
	3.00	A-TMGB-330-3	A-TMGB-380-3
	4.00	A-TMGB-330-4	A-TMGB-380-4

Temporary abutments

PEEK temporary abutments with titanium base					
	ø 3.30	ø 3.80	ø 4.25	ø 5.00	
	A-MPSC-330	-	-	-	repositionable, straight emergence
	A-MPSCR-330	A-MPSCR-380	AS-MPSCR-425	AS-MPSCR-500	repositionable, anatomical emergence

Fixation screw VM2-180 included.

SIMPLE titanium temporary abutments					
	ø 3.30	ø 3.80	ø 4.25	ø 5.00	
	A-MPSCI-330-EX	A-MPSCI-380-EX	AS-MPSCI-425-EX	AS-MPSCI-500-EX	repositionable, straight emergence
	A-MPSCI-330	A-MPSCI-380	AS-MPSCI-425	AS-MPSCI-500	rotating, straight emergence

Fixation screw VM2-180 included.

Titanium temporary abutments with castable sleeve					
	ø 3.30	ø 3.80	ø 4.25	ø 5.00	
	A-CTI-330	A-CTI-380	-	-	rotating
	A-CCI-S	A-CCI-S	-	-	spare castable PMMA sleeve, fixation screw not included

Fixation screw VM2-180 included.

Impression and model phase

Pick-up transfers					
	ø 3.30	ø 3.80	ø 4.25	ø 5.00	
	A-TRA-330	A-TRA-380	AS-TRA-425	AS-TRA-500	straight emergence
	A-TRAR-330	A-TRAR-380	AS-TRAR-425	AS-TRAR-500	anatomical emergence

Fixation screw VTTRA2-180-15 included.

Pull-up transfers					
	ø 3.30	ø 3.80	ø 4.25	ø 5.00	
	A-TRAP-330	-	-	-	straight emergence
	A-TRARP-330	A-TRARP-380	AS-TRARP-425	AS-TRARP-500	anatomical emergence

Closed-tray transfers					
	ø 3.30	ø 3.80	ø 4.25	ø 5.00	
	A-TRAS-330	-	-	-	straight emergence
	A-TRARS-330	A-TRARS-380	AS-TRARS-425	AS-TRARS-500	anatomical emergence

Fixation screw VTTRA2-180-10 included.

Analog		
	A-ANA-330	ø 3.30
	A-ANA-380	ø 3.80
	SH-ANA-425	ø 4.25
	SH-ANA-500	ø 5.00

Abutments

Premade abutments

Premade abutments for direct screwing					
	transgingival h	ø 3.30 ●	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
		A-MDAD-330	A-MDAD-380	-	-

Premade straight abutments					
	transgingival h	ø 3.30 ●	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
	1.00	A-MD-330-1	A-MD-380-1	AS-MD-425-1	AS-MD-500-1
	2.00	A-MD-330-2	A-MD-380-2	AS-MD-425-2	AS-MD-500-2
	4.00	A-MD-330-4	A-MD-380-4	AS-MD-425-4	AS-MD-500-4
	1.00	A-MDR-330-1	A-MDR-380-1	AS-MDR-425-1	AS-MDR-500-1
	2.00	A-MDR-330-2	A-MDR-380-2	AS-MDR-425-2	AS-MDR-500-2
	4.00	A-MDR-330-4	A-MDR-380-4	AS-MDR-425-4	AS-MDR-500-4

Fixation screw VM2-180 included.

Premade angled abutments					
	transgingival h	ø 3.30 ●	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
	1.75	A-MA15-330	-	-	-
	1.80	A-MAR15-330	A-MAR15-380	AS-MAR15-425	AS-MAR15-500
	1.80	-	A-MAR25-380	AS-MAR25-425	-

Fixation screw VM2-180 included.

Preparable abutments

Straight preparable abutments					
	ø 3.30 ●	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●	
	A-MF-330	A-MF-380	AS-MF-425	AS-MF-500	straight emergence
	A-MFR-330	A-MFR-380	AS-MFR-425	AS-MFR-500	anatomical emergence

Fixation screw VM2-180 included.

Pre-cut preparable abutments					
	ø 3.30 ●	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●	
	A-MFP-330	-	-	-	straight emergence
	A-MFPR-330	A-MFPR-380	AS-MFPR-425	AS-MFPR-500	anatomical emergence

Fixation screw VM2-180 included.

SIMPLE preparable abutments					
	ø 3.30 ●	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●	
	A-MFS-330	A-MFS-380	AS-MFS-425	AS-MFS-500	very wide emergence

Fixation screw VM2-180 included.

Interceptive technique preparable abutments		
	A-MFTI-380	ø 3.80 ●

Fixation screw VM2-180 included.

Prosthetic solutions for vertical technique

Components for vertical technique		
	unique platform	
	A-TMG-MEFL-330	Gr. 5 titanium healing abutments
	A-MEFL-330	Gr. 5 titanium preparable abutments
	CAP-MEFL-5	impression cap

Preparable abutments: fixation screw VM2-180 included.

XA prosthetic range

Abutments for XA cemented restorations					
	transgingival h	ø 3.30 ●	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
	1.00	SH-MD-F-330-1	SH-MD-F-380-1	SH-MD-F-425-1	-
	2.00	SH-MD-F-330-2	SH-MD-F-380-2	SH-MD-F-425-2	-

Fixation screw L-VMS-180 included.

Components for XA cemented restorations					
	ø 3.30 ●	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●	
	SH-ANA-MD-F-330	SH-ANA-MD-F-380	SH-ANA-MD-F-425	-	analog
	SH-TRA-MD-F-330	SH-TRA-MD-F-380	SH-TRA-MD-F-425	-	transfer
	SH-CP-MD-F-330	SH-CP-MD-F-380	SH-CP-MD-F-425	-	PEEK cap

Abutments for XA screw retained restorations			
	transgingival h	ø 3.30 ●	ø 3.80 ●
	4.50	A-ABU-F-TS-330-4	SH-ABU-F-TS-380-4
	5.50	A-ABU-F-TS-330-5	SH-ABU-F-TS-380-5
	6.50	A-ABU-F-TS-330-6	SH-ABU-F-TS-380-6

Components for XA screw retained restorations		
	unique platform	
	SH-ANABU-F-380	analog
	SH-TRABU-F-380	transfer
	SH-CCABU-F-380	castable sleeve, repositionable
	SH-CCABU-F-380-ROT	castable sleeve, rotating
	SH-CG-ABU-F-380	PEEK healing cap
	SH-CTABU-F-380	Gr. 5 titanium sleeve, repositionable
	SH-CTABU-F-380-ROT	Gr. 5 titanium sleeve, rotating
	CAP-TS-DEF	Conoweld cap for luting

Transfer: fixation screw SH-VTRABU-F-200 included. Castable and titanium sleeves: fixation screw A-PLAIN-VP200 included.

Castable abutments with a metal base

Castable abutments with a metal base					
	ø 3.30 ●	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●	
	A-UCR-330-EX	A-UCR-380-EX	AS-UCR-425-EX	AS-UCR-500-EX	gold alloy "1" base, repositionable
	A-UCTR-330-EX	A-UCTR-380-EX	-	-	Gr. 5 titanium base, repositionable
	A-UCRCO-330-EX	A-UCRCO-380-EX	AS-UCRCO-425-EX	AS-UCRCO-500-EX	cobalt chrome base, repositionable
	A-UCR-330	A-UCR-380	AS-UCR-425	AS-UCR-500	gold alloy "1" base, rotating
	A-UCRCO-330	A-UCRCO-380	AS-UCRCO-425	AS-UCRCO-500	cobalt chrome base, rotating
	A-CCUCR-330	A-CCUCR-380	AS-CCUCR-425	AS-CCUCR-500	spare castable PMMA sleeve, screw not included

Fixation screw VM2-180 included.

Abutments

P.A.D. abutments (Disparallel Screw retained Prosthesis)

Straight P.A.D. abutments					
	transgingival h	ø 3.30	ø 3.80	ø 4.25	ø 5.00
	1.50	A-PAD-AD330-15	A-PAD-AD380-15	AS-PAD-AD425-15	AS-PAD-AD500-15
	3.00	A-PAD-AD330-30	A-PAD-AD380-30	AS-PAD-AD425-30	AS-PAD-AD500-30
	4.00	A-PAD-AD330-40	A-PAD-AD380-40	AS-PAD-AD425-40	AS-PAD-AD500-40

Angled P.A.D. abutments						
	transgingival h	ø 3.30 	ø 3.80 	ø 4.25 	ø 5.00 	
	3.00	A-PAD-AA330-173	A-PAD-AA380-173	AS-PAD-AA425-173	AS-PAD-AA500-173	17°
	5.00	A-PAD-AA330-175	A-PAD-AA380-175	AS-PAD-AA425-175	AS-PAD-AA500-175	17°
	3.00	A-PAD-AA330-303	A-PAD-AA380-303	AS-PAD-AA425-303	AS-PAD-AA500-303	30°
	5.00	A-PAD-AA330-305	A-PAD-AA380-305	AS-PAD-AA425-305	AS-PAD-AA500-305	30°

PAD-VM-180 fixation screw included.

P.A.D. components for overstructures

Protection caps		
	unique platform	
	PAD-CG	in Gr. 5 titanium
	PAD-CGP	in PEEK

Fixation screw PAD-VTRAL-140 included.

Impression and model phase		
	unique platform	
	PAD-CAP	POM caps, rotating
	PAD-CAP-EX	POM caps, repositionable
	PAD-TRA	Gr. 5 titanium Pick-up transfer, rotating
	PAD-TRA-EX	Gr. 5 titanium Pick-up transfer, repositionable
	PAD-ANA	Gr. 5 titanium P.A.D. abutment analog

Pick-up transfer: fixation screw PAD-VTRAL-140 included.

Castable abutments with a metal base		
	unique platform	
	PAD-UC	PMMA castable sleeve, gold alloy "1" base rotating
	PAD-UCRCO	PMMA castable sleeve, cobalt chrome base rotating

Fixation screw PAD-VP-140 included.

P.A.D. components for prosthesis relining and luting technique		
	unique platform	
	PAD-CP	PEEK sleeve for relining, rotating
	PAD-CP-EX	PEEK sleeve for relining, repositionable
	PAD-CT	Gr. 5 titanium sleeve for relining and luting rotating
	PAD-CT-EX	Gr. 5 titanium sleeve for relining and luting repositionable
	PAD-CCEM	PMMA castable sleeve for Gr. 5 titanium sleeves, fixation screw not included
	PAD-CC	PMMA castable sleeve, rotating
	PAD-CC-EX	PMMA castable sleeve, repositionable

PEEK, Gr.5 titanium and PMMA: fixation screw PAD-VP-140 included.

P.A.D. prostheses for D.P.F. technique		
	unique platform	
	PAD-CT-LV	spare titanium sleeve
	PAD-CC-LV	spare castable centring device
	PAD-TR-LV	spare anti-escape plug
	PAD-ORING-LV	spare protection O-ring
	PAD-LV	complete pack of the prosthetic components for the "D.P.F." technique
	BARC	castable bar, l 5 cm, ø 2.20 mm

Titanium sleeve: fixation screw PAD-VP-140 included.

Intermediate abutments

Intermediate abutments					
	transgingival h	ø 3.30 ●	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
	1.00	A-ABU-330-1	A-ABU-380-1	-	-
	2.00	A-ABU-330-2	A-ABU-380-2	-	-

Fixation screw A-VABU-180 included.

Components for intermediate abutments				
	ø 3.30 ●	ø 3.80 ●		
	A-TRABU-330	A-TRABU-380	Pick-up transfer	
	A-ANABU-330	A-ANABU-380	analog	
	-	A-CTABU-380-ROT	Gr. 5 titanium sleeve	
	A-CCABU-330-ROT	A-CCABU-380-ROT	PMMA spare castable sleeve, fixation screw not included	

Pick-up transfer: A-VTRABU-180 screw included. Titanium sleeve: A-VABU-180 screw included.

PLAIN abutments

PLAIN abutments					
	transgingival h	ø 3.30 ●	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
	2.00	A-PLAIN-ABU330-2	A-PLAIN-ABU380-2	AS-PLAIN-ABU425-2	AS-PLAIN-ABU500-2
	3.00	A-PLAIN-ABU330-3	A-PLAIN-ABU380-3	AS-PLAIN-ABU425-3	AS-PLAIN-ABU500-3
	4.00	A-PLAIN-ABU330-4	A-PLAIN-ABU380-4	AS-PLAIN-ABU425-4	AS-PLAIN-ABU500-4

Components for PLAIN abutments					
	ø 3.30 ●	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●	
	A-PLAIN-CG330	A-PLAIN-CG380	A-PLAIN-CG425	A-PLAIN-CG500	healing cap
	A-PLAIN-CC330	A-PLAIN-CC380	A-PLAIN-CC425	A-PLAIN-CC500	castable sleeve
	A-PLAIN-CT330	A-PLAIN-CT380	A-PLAIN-CT425	A-PLAIN-CT500	Gr. 5 titanium sleeve
	A-PLAIN-TRA-330	A-PLAIN-TRA-380	A-PLAIN-TRA-425	A-PLAIN-TRA-500	Pick-up transfer
	A-PLAIN-ANA-330	A-PLAIN-ANA-380	A-PLAIN-ANA-425	A-PLAIN-ANA-500	analog

Castable and titanium sleeve: fixation screw A-PLAIN-VP200 included. Pick-up transfer: fixation screw A-PLAIN-VTRA200 included.

CAD-CAM prosthetic solutions

Scanbodies

Laboratory scanbodies in Ergal			
	ø 3.30 ●	ø 3.80 ●	
A-CAMETRA330		A-CAMETRA380	to transfer implant connection position
unique platform			
			to transfer P.A.D. abutments connection position
PAD-CAMETRA500			

Scanbody for implants: fixation screw VM2-180 included. Scanbody for abutments: fixation screw PAD-VP-140 included.

Scanbodies for intraoral use in PEEK with titanium base			
	ø 3.30 ●	ø 3.80 ●	
A-INT-CAMTRA330		A-INT-CAMTRA380	to transfer implant connection position
unique platform			
			to transfer P.A.D. abutments connection position
PAD-INT-CAMTRA			

Scanbody for implants: fixation screw VM2-180 included. Scanbody for abutments: fixation screw PAD-VP-140 included.

Scanbodies for laboratory and intraoral use in titanium			
	ø 3.30 ●	ø 3.80 ●	
A-INT-CAMTRA330-L		A-INT-CAMTRA380-L	to transfer implant connection position
unique platform			
			for P.A.D. abutment
PAD-INT-CAMTRA-L			
			for PLAIN abutment
A-INT-CAMTRA-PLAIN-L			
			for XA abutment
XA-INT-CAMTRA-L			

Scanbody for implants: fixation screw VM2-180 included. Scanbody for P.A.D. abutments: fixation screw PAD-VP-140 included. Scanbody for PLAIN and XA abutments: fixation screw A-PLAIN-VP200 included.

Titanium bases for luting technique

Standard T-Connect						
	cone height	ø 3.30 ●	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●	
	4.00	A-BASTZR-S-330-4	A-BASTZR-S-380-4	AS-BASTZR-S-425-4	AS-BASTZR-S-500-4	repositionable
	6.00	A-BASTZR-S-330-6	A-BASTZR-S-380-6	AS-BASTZR-S-425-6	AS-BASTZR-S-500-6	
	4.00	A-BASTZR-M-330-4	A-BASTZR-M-380-4	AS-BASTZR-M-425-4	AS-BASTZR-M-500-4	rotating
	6.00	A-BASTZR-M-330-6	A-BASTZR-M-380-6	AS-BASTZR-M-425-6	AS-BASTZR-M-500-6	

Fixation screw VM2-180 included.

T-Connect with 4.00 mm cone height and different transgingival heights						
	transgingival h	ø 3.30 ●	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●	
	1.00	A-TCR-S-330-4-1	A-TCR-S-380-4-1	-	-	repositionable
	2.00	A-TCR-S-330-4-2	A-TCR-S-380-4-2	-	-	
	3.00	A-TCR-S-330-4-3	A-TCR-S-380-4-3	-	-	
	1.00	A-TCR-M-330-4-1	A-TCR-M-380-4-1	-	-	rotating
	2.00	A-TCR-M-330-4-2	A-TCR-M-380-4-2	-	-	
	3.00	A-TCR-M-330-4-3	A-TCR-M-380-4-3	-	-	

Fixation screw VM2-180 included.

Sleeves for wax-up modeling					
	cone height	ø 3.30 ●	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
	4.00	A-CCBAS-330-4	A-CCBAS-380-4	A-CCBAS-425-4	A-CCBAS-500-4
	6.00	A-CCBAS-330-6	A-CCBAS-380-6	A-CCBAS-425-6	A-CCBAS-500-6

Titanium sleeves for P.A.D. abutments

PAD-BASTZR

Fixation screw PAD-VP-140 included.

Titanium bases for luting chairside technique

Chairside T-Connect supports with 5.00 mm h cementation cone					
	transgingival height	ø 3.30 ●	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
	0.50	A-TCS-330S-C	A-TCS-380S-C	AS-TCS-425S-C	AS-TCS-500S-C
	1.50	A-TCS-330S-C-15	A-TCS-380S-C-15	-	-
	2.50	A-TCS-330S-C-25	A-TCS-380S-C-25	-	-

Fixation screw VM2-180 included.

For further information please consult the prosthetic manual MP-IMP-PKS-ONE-E.

Conometric prosthesis

Straight conometric abutments

transgingival h	ø 3.30	ø 3.80
0.50	A-MD-TS-330-05	-
1.00	A-MD-TS-330-10	A-MD-TS-380-10
2.00	A-MD-TS-330-20	A-MD-TS-380-20
3.00	A-MD-TS-330-30	A-MD-TS-380-30

Fixation screw VM2-180 included.

Components for conometric abutments

unique platform	
ANA2-MD-TS	analog
CAP2-TS-DEF	cap for partial or full-arch fixed prosthesis
CAP2-TS-IND	non-rotating cap for single element
CAP2-TS-IMP	cap for impression taking
CAP2-TS-REM	cap for the removal
CAP2-TS-PIN	cap with pin

Angled conometric abutments

ø 3.30	ø 3.80	angulation
A-MA07-TS-330-2	A-MA07-TS-380-2	7.5°
A-MA15-TS-330-2	A-MA15-TS-380-2	15°

Fixation screw VM2-180 included.

Clinical and laboratory prosthetic screws

Standard fixation screws

VM2-180	
VM2-180-10	10 pcs packaging

Conical support screws

L-VMS-180	
L-VMS-180-10	10 pcs packaging

Screws for Pick-up transfer

VTRA2-180-15	h 15.00
VTRA2-180-20	h 20.00
VTRA2-180-MAN	for manual screwing

Screws for closed tray transfer

VTRA2-180-10	h 10.00
--------------	---------

Recommended torque for tightening prosthetic components

8-10 Ncm	surgical cover screws, healing abutments
20-25 Ncm	prosthetic screws
25-30 Ncm	direct screwing on implants prosthetic components
8-10 Ncm	transfer screws

Screws for P.A.D. abutments

PAD-VM-180	for P.A.D. abutments
PAD-VM-180-10	10 pcs packaging
PAD-VTRAL-140	for Pick-up transfer, long
PAD-VTRA-140	for Pick-up transfer, standard
PAD-VTRAL-140-MAN	for Pick-up transfer, manual screwing
PAD-VP-140	for prosthetic components on P.A.D. abutments
PAD-VP-140-10	10 pcs packaging

Screws for intermediate abutments

A-VABU-180	for intermediate abutments
A-VTRABU-180	for intermediate abutments Pick-up transfer

Screws for PLAIN abutments

A-PLAIN-VP200	for PLAIN abutments
A-PLAIN-VP200-10	10 pcs packaging
A-PLAIN-VTRA200	for PLAIN abutments Pick-up transfer

Screws for CAD CAM customized prosthesis

A-CAMTVABU180	for zirconium oxide single crowns and ECHO overstructures
PAD-VCAM140	for zirconium oxide ECHO overstructures on P.A.D. abutments
CAMPRON205-10	spare long ring washers for the head of the screw

Warning

Fixation screws are always included in the packaging of each prosthetic component, unless otherwise indicated.

Please use laboratory screws for laboratory phases and keep the screw included with the prosthetic component for the final fastening in the patient's mouth.

Solutions for overdenture

Locator abutments and accessories

Locator abutments		
	transgingival height	unique platform
	1.00	1670
	2.00	1671
	3.00	1672
	5.00	1674
Accessories for Locator abutments		
	unique platform	
	8530	4 aluminium analogs
	8505	4 aluminium transfers 4 retainers cod. 8515 included
	8517	4 parallelism pins in nylon
	9530	steel plate AISI 316 L
Locator Core Tool		
	8393	complete Locator Core Tool
	8397	steel push rod for inserting the retainers in the caps
	8390	steel driver for abutment screwing/ unscrewing
	8394	retention jacket for the driver, to transfer the Locator abutments into the oral cavity

Ball attachments and accessories

Ball attachments		
	transgingival height	unique platform
	1.00	A-AS-330-1
	2.00	A-AS-330-2
	4.00	A-AS-330-4
Caps for ball attachments		
	unique platform	
	CAP-TFL-1	polyamide cap
	CONT-CAP-TFL-1	steel container for polyamide cap with outer ø 4.80 mm
	CAP-TIT-1	Gr. 5 titanium cap complete with cap in two parts, titanium retention spring and plastic mounting ring
	AN-CAP-TIT-1	spare plastic ring
	MOL1-CAP-TIT-1	spare retention spring for titanium caps, average hardness
	MOL2-CAP-TIT-1	spare retention spring for titanium caps, soft
	CAP-1	gold alloy "2" cap, complete with 3 plastic positioning rings and a tin spacer for the laboratory

Retainers and spacer rings for Locator abutments		
	unique platform	
	8515	low retention capacity for impression taking, black, in polyethylene, 4 pcs packaging
	8524	retention 5 lb = 2268 g, in nylon, disparallelisms up to 10° per side, 4 pcs packaging
	8527	retention 3 lb = 1361 g, in nylon, disparallelisms up to 10° per side, 4 pcs packaging
	8529	retention 1,5 lb = 680 g, in nylon, disparallelisms up to 10° per side, 4 pcs packaging
	8547	retention 4 lb = 1814 g, in nylon, disparallelisms up to 20° per side, 4 pcs packaging
	8915	retention 2 lb = 907 g, in nylon, disparallelisms up to 20° per side, 4 pcs packaging
	8548	retention 1 lb = 450 g, in nylon, disparallelisms up to 20° per side, 4 pcs packaging
	8514	spacer rings in silicon rubber, 20 pcs packaging
	8519-2	kit: 2 Gr. 5 titanium caps, 2 spacer rings, 2 black retainers and 2 nylon retainers for each of the 3 different retention capacities
	8540-2	kit: 2 Gr. 5 titanium caps, 2 spacer rings, 2 black retainers and 2 nylon retainers for each of the 3 different retention capacities meant for disparallelisms
	8550-2	kit: 2 Gr. 5 steel caps, 2 spacer rings, 2 black retainers and 2 nylon retainers for each of the 3 different retention capacities

Analog for ball attachments	
	unique platform
	ANAS



rev. 03-22



 **Sweden & Martina S.p.A.**
Via Veneto, 10
35020 Due Carrare (PD), Italy

Authorized Distributor
Sweden & Martina Inc.
4700 S Mill Ave, Ste B-16
PO Box 23748 - Tempe, AZ 85282
Toll free 1-844-8MARTINA (1-844-862-7846)
info.us@sweden-martina.com

www.sweden-martinainc.com

The implants, standard prosthetic components and surgical instruments contained in this catalogue are Medical devices and are manufactured by Sweden & Martina S.p.A. They are conform to the ISO 9001 and ISO 13485 standards and are certified with the CE Mark (Class I), CE 0051 (ClassIr) and CE 0476 mark (Class IIA and class IIB) in compliance with Regulation (EU) Medical Devices n.2017/745. They are conform to the QSR and 21 CFR part 820 and are cleared by FDA.

The Sweden & Martina plant manufactures Medical Devices in compliance with the CGMPs in force in the USA and in other countries worldwide.



Some products may not be regulatory/released for sale in all markets.
All trademarks herein are the property of Sweden & Martina S.p.A. unless otherwise indicated.
This material is intended for laboratories and clinicians and is not intended for patient distribution.
This material is not to be redistributed, duplicated, or disclosed without the express written consent of Sweden & Martina S.p.A.
For additional product information, including indications, contraindications, warnings, precautions, and potential adverse effects, see Sweden & Martina S.p.A. website.

The contents are updated at the time of publication. Check with the company for any subsequent updates.