


PRAMAWPOWER

Product Overview





Sweden & Martina



Dental Arena











PRAMAW POWER

Table of contents

4	PramaW Slim Power implants
5	PramaW Power implants
6	PramaW RF Power implants
7	PramaW RF SL Power implants
8	Ergonomics of surgical procedures
10	PramaW Power Surgical Kit
12	Surgical Kits for Guided Surgery
16	Guided surgery accessories
18	Optional integration sets
20	Optional surgical instruments
24	Prosthetic solutions of excellence
26	Impression and model phase
27	Healing abutments
28	Prama IN healing abutments
29	Provisional posts
30	Premade posts
31	Preparable posts and Castable posts
32	P.A.D.r (Disparallel Screwed Prosthesis) prosthetic range
34	PLAIN B-Space prosthetic range
35	Conico conometric prosthesis
36	Components for digital prosthesis
40	Solutions seating on implant platform
42	Premilled
43	Overdenture solutions
46	Clinical and laboratory prosthetic screws
47	Bibliography














PramaW Slim Power implants



	SHORT NECK h 1.80 mm	REGULAR NECK h 2.80 mm
h fixture	ø 3.30 ●	
6.00 mm	-	 LAW-ZT-330-060
8.50 mm	 LASW-ZT-330-085	 LAW-ZT-330-085
10.00 mm	 LASW-ZT-330-100	 LAW-ZT-330-100
11.50 mm	 LASW-ZT-330-115	 LAW-ZT-330-115
13.00 mm	 LASW-ZT-330-130	 LAW-ZT-330-130
15.00 mm	-	 LAW-ZT-330-150
surgical fixation screw (included)	 L-VT-340	

PramaW Power implants



PRAMAW POWER			
h fixture	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
6.00 mm	 ø 3.80 ø 2.75 6.00 LAW-ZT-380-060	 ø 4.25 ø 2.90 6.00 LAW-ZT-425-060	 ø 5.00 ø 3.30 6.00 LAW-ZT-500-060
8.50 mm	 ø 3.80 ø 2.75 8.50 LAW-ZT-380-085	 ø 4.25 ø 2.90 8.50 LAW-ZT-425-085	 ø 5.00 ø 3.30 8.50 LAW-ZT-500-085
10.00 mm	 ø 3.80 ø 2.75 10.00 LAW-ZT-380-100	 ø 4.25 ø 2.90 10.00 LAW-ZT-425-100	 ø 5.00 ø 3.30 10.00 LAW-ZT-500-100
11.50 mm	 ø 3.80 ø 2.75 11.50 LAW-ZT-380-115	 ø 4.25 ø 2.90 11.50 LAW-ZT-425-115	 ø 5.00 ø 3.30 11.50 LAW-ZT-500-115
13.00 mm	 ø 3.80 ø 2.75 13.00 LAW-ZT-380-130	 ø 4.25 ø 2.90 13.00 LAW-ZT-425-130	 ø 5.00 ø 3.30 13.00 LAW-ZT-500-130
15.00 mm	 ø 3.80 ø 2.75 15.00 LAW-ZT-380-150	 ø 4.25 ø 2.90 15.00 LAW-ZT-425-150	-
surgical cover screw (included)	 L-VT-340		

PramaW RF Power implants



PRAMAW RF POWER			
h fixture	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
6.00 mm	-	 LSW-ZT-425-060	 LSW-ZT-500-060
8.50 mm	 LSW-ZT-380-085	 LSW-ZT-425-085	 LSW-ZT-500-085
10.00 mm	 LSW-ZT-380-100	 LSW-ZT-425-100	 LSW-ZT-500-100
11.50 mm	 LSW-ZT-380-115	 LSW-ZT-425-115	 LSW-ZT-500-115
13.00 mm	 LSW-ZT-380-130	 LSW-ZT-425-130	 LSW-ZT-500-130
15.00 mm	 LSW-ZT-380-150	 LSW-ZT-425-150	-
surgical cover screw (included)	 L-VT-340		


PramaW RF SL Power implants



PRAMAW RF SL POWER			
h fixture	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
6.00 mm	-	-	-
8.50 mm	 LSW-ZT-380SL-085	 LSW-ZT-425SL-085	 LSW-ZT-500SL-085
10.00 mm	 LSW-ZT-380SL-100	 LSW-ZT-425SL-100	 LSW-ZT-500SL-100
11.50 mm	 LSW-ZT-380SL-115	 LSW-ZT-425SL-115	 LSW-ZT-500SL-115
13.00 mm	 LSW-ZT-380SL-130	 LSW-ZT-425SL-130	 LSW-ZT-500SL-130
15.00 mm	 LSW-ZT-380SL-150	 LSW-ZT-425SL-150	-
surgical cover screw (included)	 L-VT-340		

Ergonomics of Surgical Procedures

The **SURGICAL KITS** are comprehensive, **compact, ergonomic**, and contain all the necessary tools for the surgery of **PramaW Power implants**.



The color codes, also present on the surgical trays, guide the correct sequence of instrument usage.



Both the drills and the stops are characterized by a color code that guides operators in tool selection.



The **SURGICAL TRAYS** are made of **Radel**, a polymer that maintains its characteristics even after numerous cycles of washing and autoclave sterilization.

CYLINDRICAL or CONICAL DRILLS, based on different preparations, are supplied complete with depth stops.

The drills are made of **stainless steel** for surgical use and maintain their characteristics unchanged over time and through multiple uses.



The drills are characterized by the presence of clearly visible depth notches.

The drills are **highly sharp, centered, and perform well**. The morphology of the conical drills allows for the collection of valuable bone for regeneration needs.





The kits include a convenient **RATCHET**, which can be used both as a fixed wrench and as a dynamometric wrench for controlling the implant insertion torque and tightening of prosthetic screws.

- The torque can be adjusted from 10 to 70 Ncm;
- The torque adjustment is highly accurate because each individual ratchet is individually calibrated.

Included is a wide range of **PROSTHETIC DRIVERS** for digital use, for the dynamometric ratchet and for micromotor.

The range is suitable for different oral access needs, from extra-short tools to longer instruments for Toronto-type prosthetic rehabilitations.

The slightly tapered coupling between the driver and the screw allows for adequate retention, ensuring secure transport of the screw in the oral cavity.




The surgical procedure involves a mountless approach that ensures visibility and ergonomics with the use of the patented **EASY INSERT** drivers:

- it can be easily removed from implant wells even in the presence of high insertion torque;
- the dodecagonal design prevents deformation of the connection edges, ensuring the maintenance of precision and prosthetic stability;
- available in the kit in various configurations, both for use with a micromotor and the ratchet;
- the instruments have a hexagon on the shaft, aligned with the connection hexagon, to facilitate rotational orientation of the fixture.

PramaW Power surgical kit

Complete surgical kit containing all the necessary tools for the insertion of PramaW Power implants (code ZPRAMA-INT).
All instruments can be purchased individually as spare parts.

- 1 Precision drill**
FS-230
- 2 Parallelism pins**
PPTL-2-28
- 3 Intermediate drill**
• FGS-340/425
- 4 Pilot drill***
• FPT3-200-LXS
- 5 Stops for pilot drills**
• STOP4-200-085
• STOP4-200-100
• STOP4-200-115
• STOP4-200-130
• STOP4-200-150
- 6 Intermediate cylindrical drill**
ø 2.50 mm*
• FFT3-250-LXS
- 7 Stop for intermediate cylindrical drill**
ø 2.50 mm
• STOP4-250-085
• STOP4-250-100
• STOP4-250-115
• STOP4-250-130
• STOP3-250-150
- 8 Intermediate conical drill**
• SE-FK250
- 9 Stops for Intermediate conical drill**
CSR-STOP-2028-085
CSR-STOP-2028-100
CSR-STOP-2028-115
CSR-STOP-2028-130
CSR-STOP-2028-150
- 10 Final cylindrical drill***
• FFT3-280-LXS for implants ø 3.30 mm
• FFT3-300-LXS for implants ø 3.80 mm
• FFT3-340-LXS for implants ø 4.25 mm
• FFT3-425-LXS for implants ø 5.00 mm
- 11 Stops for conical drills**
FFT3-280-LXS
• STOP4-280-085
• STOP4-280-100
• STOP4-280-115
• STOP4-280-130
• STOP4-280-150
- 12 Stops for cylindrical drills**
FFT3-300-LXS
• STOP4-300-085
• STOP4-300-100
• STOP4-300-115
• STOP4-300-130
• STOP4-300-150
- 13 Stops for cylindrical drills**
FFT3-340-LXS
• STOP4-340-085
• STOP4-340-100
• STOP4-340-115
• STOP4-340-130
• STOP4-340-150
- 14 Stops for cylindrical drills**
FFT3-425-LXS
• STOP4-425-085
• STOP4-425-100
• STOP4-425-115
• STOP4-425-130
• STOP4-425-150
- 15 Conical drills**
• SH-FK380-085
• SH-FK380-100
• SH-FK380-115
• SH-FK380-130
• SH-FK380-150
- 16 Replies**
• SH-380-085-RP
• SH-380-100-RP
• SH-380-115-RP
• SH-380-130-RP
• SH-380-150-RP
- 17 Conical drills**
• SH-FK425-085
• SH-FK425-100
• SH-FK425-115
• SH-FK425-130
• SH-FK425-150
- 18 Replies**
• SH-425-085-RP
• SH-425-100-RP
• SH-425-115-RP
• SH-425-130-RP
• SH-425-150-RP
- 19 Conical drills**
• SH-FK500-085
• SH-FK500-100
• SH-FK500-115
• SH-FK500-130
• SH-FK500-150
- 20 Replies**
• SH-500-085-RP
• SH-500-100-RP
• SH-500-115-RP
• SH-500-130-RP
• SH-500-150-RP



The Prama Power implants can also be placed with the following kit: ZPREMIUM-ONE-INT.

Prama Power RF and Prama Power RFSL implants can also be placed with the following kit: ZSHELTAP-INT.

PSlim implants are compatible with all the previously mentioned and available surgical kits.

L-TRAYL-INT instrument carrying case in Radel.

For surgical procedures and related warnings, it is recommended to refer to the surgical manual MC-IMP-PRAMA-E (downloadable from the website www.sweden-martina.com or by scanning the QR CODE next to it).



- 21 Stops for conical drills**
- SH-STOP4-FK380
 - SH-STOP4-FK425
 - SH-STOP4-FK500

- 22 Countersink drills**
- A-FCC-380
 - A-FCC-425
 - A-FCC-500

- 23 Easy Insert Drivers for implant insertion**
- EASYC4-EX230-CA

EASYL4-EX230-CA

EASY4-EX230-EX

- 25 Torque control ratchet**
- CR15-KIT

- 24 Digital screwdriver**
- HSMXS-20-DG HSM-20-DG HSML-20-DG

- Prosthetic screwdriver**
- HSM-20-CA HSM-20-EX HSML-20-EX

- 26 Extensions and adapter**
- PROF-CAL3
Extension for surgical drills

BPM-15
Extensions for instruments with hexagonal connector

B-AVV-CA3
Mechanical adapter for instruments with hexagonal connector


- 27 Manual adapter for instruments with hexagonal connector**
- AVV3-MAN-DG

- 28 Adapter for instruments with hexagonal connector**
- AVV-CA-DG-EX

- 29 Conical bone taps**
- SH-MS-380-CA
 - SH-MS-425-CA
 - SH-MS-500-CA


- 30 Cylindrical bone taps**
- SM-MS-330
 - A-MS-380
 - A-MS-425
 - A-MS-500

- 31 Pins for neck positioning**
- L-PP-085
 - L-PP-100
 - L-PP-115
 - L-PP-130
 - L-PP-150



X-ray templates included in the kit, they can be reordered separately.


PramaW Power, PramaW RF Power










L-L100	real size
L-L120	20% larger size
L-L130	30% larger size

Surgical kit for guided surgery with PramaW Power implants

Complete surgical kit containing all the necessary instruments for guided surgery with PramaW Power implants (code ZGS-PRO-A-INT). In the tray there are a number of instrument slots referring to different systematics. These instruments are therefore not included in the standard kit.



Optional drills for the preparation of implants h 18.00 mm.



For surgical procedures and related warnings, it is recommended to refer to the surgical manual MC-UD-EPLANPROA-E (downloadable from the website www.sweden-martina.com or by scanning the QR CODE next to it).



Optional drills for preparing 6.00 mm high implants can be ordered also cumulatively as a supplementary set with the code GS-PROA-INTEGRA-060.

5

Optional drills for h 7.00 mm
ø 3.80 mm implants
GS-F200-060-415
GS-F250-060-415
GS-F280-060-415
GS-F300-060-415



6

Optional drills for h 7.00 mm
ø 4.25 and 5.00 mm implants
GS-F200-060-550
GS-F300-060-550
GS-F340-060-550
GS-F425-060-550



7

Intermediate drills
for ø 4.15 mm sleeve
GS-F200-085-415
GS-F200-100-415
GS-F200-115-415
GS-F200-130-415
GS-F200-150-415



GS-F250-085-415
GS-F250-100-415
GS-F250-115-415
GS-F250-130-415
GS-F250-150-415

8

Final drills for
ø 4.15 mm sleeve
GS-F280-085-415
GS-F280-100-415
GS-F280-115-415
GS-F280-130-415
GS-F280-150-415




GS-F300-085-415
GS-F300-100-415
GS-F300-115-415
GS-F300-130-415
GS-F300-150-415

9

Intermediate drills for
ø 5.50 mm sleeve
GS-F200-085-550
GS-F200-100-550
GS-F200-115-550
GS-F200-130-550
GS-F200-150-550



GS-F300-085-550
GS-F300-100-550
GS-F300-115-550
GS-F300-130-550
GS-F300-150-550



11

Bone taps
GS-MS-330 per impianto ø 3.30 mm
GS-MS-A380 per impianto ø 3.80 mm
GS-MS-A425 per impianto ø 4.25 mm
GS-MS-A500 per impianto ø 5.00 mm



10

Final drills for ø 4.25 mm implants
with ø 5.50 mm sleeve
GS-F340-085-550
GS-F340-100-550
GS-F340-115-550
GS-F340-130-550
GS-F340-150-550



Final drills for ø 5.00 mm implants
with ø 5.50 mm sleeve
GS-F425-085-550
GS-F425-100-550
GS-F425-115-550
GS-F425-130-550
GS-F425-150-550

Surgical kit for guided surgery with PramaW RF and PramaW RF SL Power implants

Complete surgical kit containing all the necessary instruments for guided surgery with Prama RF and Prama RF SL Power implants (code ZGS-PRO-S-INT). In the tray there are a number of instrument slots referring to different systematics. These instruments are therefore not included in the standard kit.

- 1** Mucotomes
GS-MUC-415
GS-MUC-550



- 2** Bone profilers
GS-LC-415
GS-LC-550



- 3** Initial drills and progressive drill
GS-F200-415
GS-F200-550



Optional drills for the preparation of implants h 6.00 mm, sleeve 415:

- 18** GS-F200-060-415
GS-FK410U-060-415
GS-FK410-060-415



Optional drills for the preparation of implants h 6.00 mm, sleeve 550:

- 17** GS-F200-060-550
GS-FK410U-060-550
GS-FK500U-060-550
GS-FK500-060-550



Can be ordered also cumulatively as a supplementary set with the code GS-PROS-INTEGRA-060.


- 16** Adapter for instruments with hexagonal connector
AVV3-MAN-DG



- 15** Torque control ratchet
CRI5-KIT




- 14** Pins and drill for pins
GS-FPN-148 GS-PIN (4 pcs)



- 13** Drivers for mounters
EASY4-EX250-EX EASYC4-EX250-CA




- 12** Countersink drills
GS-FCS-SH380
GS-FCS-SH425
GS-FCS-SH500
GS-FCS-SE380
GS-FCS-SE425
GS-FCS-SE500





15

Guided surgery accessories





implant ø				
3.30 ●		3.80 ●		4.25 ● - 5.00 ●
sleeve		 5.00 GS-B415	 5.00 GS-B415	 5.00 GS-B550
		 4.00 GS-B415-EX-6*	 4.00 GS-B415-EX-6*	 4.00 GS-B550-EX-6*
mounter	SHORT NECK	-	 7.20 GS-MOU2-L415-052	 7.20 GS-MOU2-L500-052
	REGULAR NECK	 9.00 GS-MOU-A330	 6.20 GS-MOU-L415	 6.20 GS-MOU-L415  9.00 GS-MOU-L550 GS-MOU-A380SP**
	LONG NECK	-	 5.20 GS-MOU2-L415-072	 5.20 GS-MOU2-L500-072
fixation screw included		 GS-VTMOU-180**	 GS-VTMOU-180**	 GS-VTMOU-180**
easy mounter		 9.00 GS-EASY-AS-415-EX	 9.00 GS-EASY-L-415-EX	 9.00 GS-EASY-L-550-EX
		 9.00 GS-EASY-AS-415-CA	 9.00 GS-EASY-L-415-CA	 9.00 GS-EASY-L-550-CA
easy mounter with depth indication		-	 GS-EASYPE-L-415-EX	 GS-EASYPE-L-550-EX
		-	 GS-EASYPE-L-415-CA	 GS-EASYPE-L-550-CA

*The engaging guide sleeves are available in packs of 6 pieces, not sold individually. Non-engaging guide sleeves are also available in packs of 6, which can be ordered using the codes GS-B415-6 and GS-B550-6.


**It is recommended, with this mounter, to adjust the implant site preparation to ensure the torque value does not exceed 50 Ncm.

16 ***The fixation screw is included with the mounter but can also be ordered individually as a spare.

Bone taps without stop

implant ø	3.30 ●	3.80 ●	4.25 ●	5.00 ●
PramaW Power				
	GS-MSPE-A330	GS-MSPE-A380	GS-MSPE-A425	GS-MSPE-A500
Prama RF Power, Prama RF SL Power				
		GS-MSPE-S380	GS-MSPE-S425	GS-MSPE-S500

Countersink drills without stop

implant ø	3.30 ●	3.80 ●	4.25 ●	5.00 ●
PramaW Power				
	-	GS-FCSPE-380	GS-FCSPE-425	GS-FCSPE-500

Sleeves for pins

6 item pack



Guidelines for PramaW Power
insertion with guided surgery



Optional integration sets

L-INTEGRA-060 set

The L-INTEGRA-060 drill set includes the components necessary for the insertion of the PramaW Power and Prama RF Power h 6.00 mm implants, the housings of which are provided inside the PramaW Power kit.

Initial drills PramaW Power h 6.00 mm*

pilot drill	Shorty intermediate drill
	
FPS-200	FGS-200/300
	-
STOPS4-200-060	-


* Instruments included in ZSHORTY-INT

Final drills PramaW Power h 6.00 mm*

ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
		
FFS-300	FFS-340	FFS-425
		
STOPS4-300-060	STOPS4-340-060	STOPS4-425-060

* Instruments included in ZSHORTY-INT

Parallelism posts



Final drills PramaW RF Power h 6.00 mm

ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
-		
-	SE-FK410	SE-FK500
-		
-	SE-STOP-FK410-060	SE-STOP-FK500-060

Reply for PramaW RF Power h 6.00 mm



Set L-INTEGRA set

L-INTEGRA set includes all the necessary components for the insertion of PramaW Slim Power $\varnothing$ 3.30 mm, the countersink drills and the posts for neck positioning. These components are included in the ultimate version of the PramaW Power surgical kit, but not in the previous version of the surgical kit, thus they are still available for those who have the previous version.

Posts for neck positioning

	h	
	8.5	L-PP-085
	10.00	L-PP-100
	11.50	L-PP-115
	13.00	L-PP-130
	15.00	L-PP-150

Countersink drills

$\varnothing$ 3.80 	$\varnothing$ 4.25 	$\varnothing$ 5.00 
		
A-FCC-380	A-FCC-425	A-FCC-500

Instruments for Prama Slim Power implant

final drill	stop	h	bone tap with hexagonal connector
		8.50	
		10.00	
		11.50	
		13.00	
		15.00	
FFT3-280-LXS		STOP4-280-085	A-MS-330
		STOP4-280-100	
		STOP4-280-115	
		STOP4-280-130	
		STOP4-280-150	

Optional surgical instruments (not included in PramaW Power surgical kit)

Long precision drill



FSL-230

Intermediate drills



FG-200/280XS

intermediate drill 200/280



FG-330/425XS

intermediate drill 330/425

Drills for distal sectors

ø 2.00



FFT5-200-LXS

ø 3.30



FFT5-330-LXS

ø 2.80



FFT5-280-LXS

ø 3.40



FFT5-340-LXS

ø 2.90



FFT5-290-LXS

ø 3.60



FFT5-360-LXS

ø 3.00



FFT5-300-LXS

ø 4.25



FFT5-425-LXS

ø 3.20



FFT5-320-LXS

ø 4.45



FFT5-445-LXS

Screwdrivers not included in the kit



HSMXL-20-EX

with hexagonal connector, extra long



HSMXS-20-CA

with right angle shank, extra short



L-HSM-EX

Full Head, with hexagonal connector, short



L-HSML-EX

Full Head, with hexagonal connector, long



L-HSMXL-EX

Full Head, with hexagonal connector, extra long



L-HSM-CA

Full Head, with right angle shank



AVV-CAP-TIT-1

instrument for maintaining the titanium cap for ball attachments



BASCC-EX

for ball attachments, with hexagonal connector



AVV2-ABUT

for straight P.A.D. abutments



8926-SW

for Locator abutment, short



8927-SW

for Locator abutment, long



DSPDCLH-24


for angled screw hole, length 24 mm



DSPDCLH-32


for angled screw hole, length 32 mm

Parallel pins




PP-2/28

Parallelism pin with one side of $\varnothing$ 2.00 mm and the other side of $\varnothing$ 2.80 mm



PPTS-2-28

Parallelism pin with depth notches, small version



PPS-2/3

Parallelism pin for short drills with notches at 5.00, 6.00, and 7.00 mm

Bone taps

$\varnothing$ 3.30

$\varnothing$ 3.80

$\varnothing$ 4.25

$\varnothing$ 5.00

Bone taps with hexagonal connector



SM-MSC-330

A-MSC-380

A-MSC-425

A-MSC-500



-

SH-MS-380

SH-MS-425

SH-MS-500

Bone taps with right angle shank



SM-MS-330-CA

A-MS-380-CA

A-MS-425-CA

A-MS-500-CA

Extraction tools for intraoral removing of implants



BC-EX230

short



BL-EX230

long

Extractors for atraumatic removal of implants



EXTR-IMP-S

small



EXTR-IMP-M

medium



EXTR-IMP-L

large

Bone profiler for P.A.D. abutments

unique platform



A-PAD-PS330-L

transgingival height 3.00, wide

A-PAD-PS330-S

transgingival height 5.00, narrow



A-PAD-GUI-PS-230

guide for profiler

Osteotomes



ø 0.20

E-OS-020-PP

osteotome ø 0.20
flat tip



ø 0.90

E-OS-090-PP

osteotome ø 0.90
flat tip



ø 1.60

E-OS-160-PC

osteotome ø 1.60
concave tip



ø 2.00

E-OS-200-PC

osteotome ø 2.00
concave tip



ø 2.40

E-OS-240-PC

osteotome ø 2.40
concave tip

Osteotomes with flat tip



h	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
8.50-10.00	SH-OS-380-100-PP	SH-OS-425-100-PP	SH-OS-500-100-PP
11.50	SH-OS-380-115-PP	SH-OS-425-115-PP	SH-OS-500-115-PP
13.00	SH-OS-380-130-PP	SH-OS-425-130-PP	SH-OS-500-130-PP
15.00	SH-OS-380-150-PP	SH-OS-425-150-PP	SH-OS-500-150-PP

Osteotomes with round tip



h	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
8.50-10.00	SH-OS-380-100-PR	SH-OS-425-100-PR	SH-OS-500-100-PR
11.50	SH-OS-380-115-PR	SH-OS-425-115-PR	SH-OS-500-115-PR
13.00	SH-OS-380-130-PR	SH-OS-425-130-PR	SH-OS-500-130-PR
15.00	SH-OS-380-150-PR	SH-OS-425-150-PR	SH-OS-500-150-PR

Torque wrench with control lever



TWL



TWL-AVV-EX

Adaptor for
torque wrench

Depth gauge



PROF3

Spare O-ring



ORING180-088

Service mounter for impact procedures



MOU-EX230

mounter



CM2

mounter
stop key

Soft tissue depth gauge




MP-LMRC-180

h 10 mm


Depth notches

The laser marking on the cutting edges and the use of stops allow for safe preparation and control of the instrument's insertion depth. The drills always create a hole longer than the desired implant. The oversizing (LS) is equal to the height of the tip of the drill being used.

LL: total length of the working part, including the tip.
LS: length of the overpreparation. This measurement must be calculated in addition to the length of the preparation hole.



The working depth is also identified on bone taps and osteotomes through laser-etched notches, as follows:



The working depth is also identified on bone taps without a stop and Countersink drills through laser markings, as follows:



Prosthetic solutions of excellence

Marginless solutions for luted crown or bridges


L-MD*: freedom to decide the decision on where to close the prosthetic crown margin. These solutions promote soft tissues thickening because they leave more space for connective tissue.

L-MDT*: these components have a screw channel that either allows use of the standard straight connecting screw or allows palatal or lingual reorientations of the dedicated angled screw to enable esthetic crown solution with no vestibular holes.

L-TCIN*: when single rehabilitation involves a terminal molar, the “Prama IN” prosthetic solution offers an effective option. By closing around the implant neck, it provides greater strength to the restoration, particularly suited to withstand the significant chewing forces typical of this area.

Additionally, the “Prama IN” solution can be a valid choice in cases where bone and soft tissue remodeling did not meet expectations due to individual immune response, or when space recovery is needed in areas with limited vertical height for mastication.

T-CUT: when choosing the PRAMA IN solution, the practical links designed for luting restorations offer adjustable height, with three different options available. This versatility allows for optimal positioning of the screw channel, either palatally or lingually, ensuring both functional and aesthetic advantages in the final rehabilitation. Their adjustability reduces the need for a large inventory, simplifying clinical procedures and stock management.



Solutions with predetermined base for luted crowns or bridges


T-UNICA: T-Connect Unica components are unique in design, featuring a standardized luting cone with a single height of 8 mm for all connections. This can be easily adjusted to 6 mm or 4 mm using dedicated cutting guides, offering versatility without the need for a large inventory. The emergence profile reaches a diameter of 4.40 mm, with a variable transmucosal height ranging from 0.80 mm to 4 mm, allowing adaptation to different clinical situations. Their adjustability simplifies both clinical procedures and stock management, making them a practical and efficient solution.

BASTZR*: standard shoulder T-Connects are available for straight screw channels and offered in different transgingival heights, allowing adaptation to various clinical situations.



CONICO

- Conometric prosthesis for fixed restorations without screws or cement.
- Easily removable: the clinician can detach the prosthesis in seconds for follow-ups and hygiene procedures.
- Suitable for restorations ranging from single crowns to full arches, thanks to the conometric technique.
- Ensures a biological seal, promoting healthy and stable tissues over time.
- Minimizes chair time during check-ups, thanks to the simplicity of the system.



P.A.D.r abutments for multiple prosthetics





particularly suitable in rehabilitation of the partially or totally edentulous arches and in situations of reduced horizontal and vertical spaces.

- Easy to insert;
- maximum freedom and ergonomics of the superstructures;
- pink anodization for better mimetism;
- their limited size makes them suitable for rehabilitations of partially or totally edentulous arches in reduced horizontal and vertical spaces;
- allow tilted implants protocols




Impression and model phase

Pick-up transfers




straight emergence			
ø 3.30 ●	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
			
A-TRA-330	L-TRA-380	L-TRA-425	L-TRA-500

A-VTRA-180-15 screw included.







Pull-up transfers

straight emergence
unique platform

A-TRAP-330

Analogs

ø 3.30 ●	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
			
A-ANA-330	L-ANA-380	L-ANA-425	L-ANA-500

Healing abutments

Healing abutments					
transgingival h	2 .00 mm	3 .00 mm	4 .00 mm	5 .00 mm	7 .00 mm
with straight emergence, unique platform	 A-TMG-330-2	 A-TMG-330-3	-	 A-TMG-330-5	 A-TMG-330-7
with straight emergence with coronal enlargement, unique platform	 A-TMTCR-330-2	 A-TMTCR-330-3	-	 A-TMTCR-330-5	-
with wide profile, unique platform	 A-TMGB-330-2	 A-TMGB-330-3	 A-TMGB-330-4		
Healing abutment					
unique platform					
titanium Gr. 5	 A-TMG-MEFL-330				

Prama IN healing abutments

(closing around the implant neck)


in Titanium						
transgingival h	closing around the neck h	ø 3.80	ø 4.25	ø 5.00		
2.50	0.50*	5.00 LW-TMG-380-25-05	5.40 LW-TMG-425-25-05	6.10 LW-TMG-500-25-05		
3.50	0.50*	5.00 LW-TMG-380-35-05	5.40 LW-TMG-425-35-05	6.10 LW-TMG-500-35-05		
2.50	1.50**	5.00 LW-TMG-380-25-15	5.40 LW-TMG-425-25-15	6.10 LW-TMG-500-25-15		
3.50	1.50**	5.00 LW-TMG-380-35-15	5.40 LW-TMG-425-35-15	6.10 LW-TMG-500-35-15		

*Only for Regular and Long neck. | *Only for Regular neck.

Provisional posts


Simple in PEEK and Titanium Gr. 5

unique platform



A-MPSC-330

Simple in PEEK, base in Titanium Gr. 5, repositionable



A-MPSCI-330-EX

Simple in Titanium Gr. 5, repositionable

A-MPSCI-330

Simple in Titanium Gr. 5, rotating

VM2-180 fixation screw included.

In Titanium with castable sleeve

unique platform



A-CTI-330

rotating

A-CCI-S

spare castable PMMA sleeve,
fixation screw not included

VM2-180 fixation screw included.

Posts for luting technique

unique platform



L-CT-340-ROT

Gr. 5 titanium posts, h 14.00 mm,
straight emergence

L-CCRCO-340-ROT

Cobalt chrome titanium posts,
h 14.00 mm, straight emergence



L-CTR-340-ROT

Gr. 5 titanium posts, h 14.00 mm,
anatomical emergence

L-CCRCOR-340-ROT




Cobalt chrome posts, h 14.00 mm,
anatomical emergence

L-VMS-180 screw included.

Premade posts


Straight, marginless			
unique platform	post h		
	15.00	A-MD-330-15	repositionable

L-VMS-180 screw included.

Angled			
unique platform	transgingival h		
	0.00	L-MA15-340	
	1.00	L-MA15-340-1	angled at 15°, repositionable
	2.00	L-MA15-340-2	
	1.00	L-MA25-340-1	angled at 25°, repositionable
	2.00	L-MA25-340-2	





L-VMS-180 fixation screw included, for L-MA15-340 post: L-VM-180 screw not included.

Preparable posts


Preparable posts	
unique platform	
straight emergence	 A-MF-330
wide emergence	 A-MFS-330

VM2-180 fixation screw included.

Castable posts

With metal base	
unique platform	
cobalt chrome base, repositionable	 A-UCRCO-330-EX
cobalt chrome base, rotating	 A-UCRCO-330
titanium base, rotating <i>only for titanium fusions</i>	 A-UCTR-330-EX
spare castable sleeve, fixation screw not included	 A-CCUCR-330

VM2-180 fixation screw included.

Fully castable	
unique platform	
repositionable, straight emergence	 A-CC-330-EX
rotating, straight emergence	 A-CC-330

VM2-180 fixation screw included.

P.A.D.r prosthetic range (Disparallel Screw-Retained Prosthesis)

Straight P.A.D.r abutments, for direct screwing on implants

transgingival h


unique platform

1.50 mm




A-PADR-AD330-15

3.00 mm



A-PADR-AD330-30

4.00 mm



A-PADR-AD330-40

Carrier PADR-VTRAL-140-M included: for intraoral screwing use the AVV2-ABUT screwdriver.

Angled abutments


unique platform

Angled at 17°
transgingival h
3.00 mm



A-PADR-AA330-173

Angled at 17°
transgingival h
5.00 mm




A-PADR-AA330-175

Angled at 30°
transgingival h
3.00 mm



A-PADR-AA330-303

Angled at 30°
transgingival h
5.00 mm




A-PADR-AA330-305

Angled at 45°
transgingival h
3.00 mm



A-PADR-AA330-453

Angled at 45°
transgingival h
5.00 mm



A-PADR-AA330-455


Screw PADR-VM-180 included.

Carrier PADR-VTRAL-140-M included: for intraoral screwing use the appropriate screwdriver from the HSM series.

P.A.D.r prosthetic components

Healing caps

Titanium healing cap



PADR-CG

PEEK healing cap




PADR-CGP

PADR-CG: PAD-VP-140 screw included. | PADR-CGP: PAD-VCGP-140 screw included.

P.A.D.r. components for prosthesis relining and luting technique

PEEK sleeve, non engaging




PADR-CP

Titanium sleeve, non engaging



PADR-CT

PMMA castable sleeve



PADR-CCEM

PMMA castable sleeve,
non engaging



PADR-CC

Titanium sleeve without shoulder, for
aesthetic reconstructions



PADR-BAS

Sleeves in PEEK, Titanium Gr. 5 and PMMA: screw PAD-VP-140 included.

Castable post with metal base

PMMA post with cobalt-chrome base,
rotating



PADR-UCRCO

PAD-VP-140 screw included.

Impression phase

snap-on cap



PADR-CAP-EX

transfer pick-up, rotating



PADR-TRA

analog




PADR-ANA

Pick up transfer: PAD-VTRAL-140 screw included.

P.A.D.r components for CAD-CAM prosthesis

Intraoral scanbody



PADR-INT-CAMTRA-L

Digital analog for P.A.D.r abutment



PADR-ANA-DIG

Sleeve for luting technique non engaging for abutment P.A.D.r with cone h 4



PADR-TC-M-4

Sleeve for luting technique non engaging for abutment P.A.D.r with cone h 8 and reduction cut



PADR-TC-M-8

Sleeve for luting technique non engaging for abutment P.A.D.r for angled screw hole



PADR-TCA-M*

Sleeve for luting technique repositionable for abutment P.A.D.r with cone h 4



PADR-TC-S-4

Sleeve for luting technique repositionable for abutment P.A.D.r with cone h 8 con tacche di riduzione



PADR-TC-S-8


Sleeve for luting technique repositionable for abutment P.A.D.r for angled screw hole



PADR-TCA-S*

P.A.D.r components for CAD-CAM prosthesis

Castable sleeve for PADR-TC-M and PADR-TC-M-8, non engaging



PADR-TC-CC-M-8

Engaging castable sleeve for PADR-TCA-M






PADR-TC-CC-S-8

Scanbody and sleeve for luting technique: PAD-VP-140 screw included.
Digital analog: VADIG-180 screw included.



*Non engaging sleeve for luting technique for P.A.D.r abutment for angled screw hole: screw PAD-VPA-140 not included.

PLAIN B-Space prosthetic range

Abutment	
transgingival h	unique platform
2.00 mm	 ø 3.30 — 2.00 A-PLAIN-ABU330-2
3.00 mm	 ø 3.30 — 3.00 A-PLAIN-ABU330-3
4.00 mm	 ø 3.30 — 4.00 A-PLAIN-ABU330-4

Components for overstructure	
	unique platform
healing cap	 ø 3.40 — 3.40 A-PLAIN-CG330
castable sleeve	 ø 3.30 — 10.00 A-PLAIN-CC330
Gr. 5 titanium sleeve	 ø 3.30 — 9.50 — 1.30 A-PLAIN-CT330
Pick-up transfer	 ø 3.30 — 12.00 A-PLAIN-TRA-330
analog	 ø 3.30 — 12.00 A-PLAIN-ANA-330

Castable and titanium sleeves: A-PLAIN-VP200 fixation screw included.
Pick-up transfers: A-PLAIN-VTRA200 fixation screw included.

Components for digital prosthesis	
	unique platform
Titanium scanbody for laboratory and intraoral use, for PLAIN B-Space abutment	 ø 3.30 — 13.00 A-INT-CAMTRA-PLAIN-L
Digital analog for PLAIN B-Space abutment	 ø 3.30 — 12.00 A-PLAIN-ANA-330-DIG
PLAIN B-Space sleeve for luting technique h. 6.00 mm	 ø 3.30 — 5.50 A-PLAIN-MDT-6
PLAIN B-Space sleeve for luting technique h. 10.00 mm	 ø 3.30 — 10.00 A-PLAIN-MDT-10

Scanbody: A-PLAIN-VP200 screw included.
Digital analog: VADIG-180 screw included.
PLAIN B-Space sleeves for luting technique: screw not included.
A-PLAIN-VP200 screw for straight hole or ABU-VMA-200 for angled screw hole available for separate ordering.

Conico conometric prosthesis

Straight

h	unique platform
0.50	 4.75 ø 3.50 0.50 A-MD-TS-330-05
1.00	 4.75 ø 3.50 1.00 A-MD-TS-330-10
2.00	 4.75 ø 3.50 2.00 A-MD-TS-330-20
3.00	 4.75 ø 3.50 3.00 A-MD-TS-330-30

VM2-180 fixation screw included.

Conico components

cap for partial or full-arch fixed prosthesis	 5.70 CAP2-TS-DEF
non-rotating cap for single element	 5.70 CAP2-TS-IND
cap for the removal	 6.30 CAP2-TS-REM
PEEK transfer cap for direct impression on Conico abutment	 9.10 CAP2-TS-IMP
Conico post analog	 16.50 ANA2-MD-TS
cap with pin, to be used for intraoral wel (welding bars available upon request)	 7.00 5.70 CAP2-TS-PIN

Angled

	unique platform
angled at 7.5°	 4.75 5.20 ø 3.50 1.50 A-MA07-TS-330-2
angled at 15°	 4.75 4.90 ø 3.50 1.50 A-MA15-TS-330-2
angled at 22.5°	 4.75 4.70 ø 3.50 2.00 A-MA20-TS-330-2

VM2-180 fixation screw included.

Conico Carrier

Parallelization carrier	 CAP2-TS-PAR
-------------------------	---

Accessories for Conico system

parallelometer	 PAR-PP
parallelometer insert	 PAR-INS-AS

Remember to always purchase the parallelometer with the insert relating to the implant connection on which it is to be used.


Components for Conico digital prosthesis

	unique platform
Titanium scanbody for laboratory and intraoral use, for Conico abutment	 10.00 ø 3.50 C-INT-CAMTRA-L
Digital analog for Conico abutment	 10.00 ANA2-MD-TS-DIG

VADIG-180 VM2-180 screw included.

Components for digital prosthesis




Components for impression and model phase - platform level

	ø 3.30 ●	ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
scanbody*				
	L-INT-CAMTRA-L	L-INT-CAMTRA-L	L-INT-CAMTRA-L	L-INT-CAMTRA-L
digital analog regular neck**				
	A-ANA-330-DIG	L-ANA-380-DIG	L-ANA-425-DIG	L-ANA-500-DIG

* VM2-18 screw included. | ** VADIG-180 screw included.

Components for impression and model phase - abutment level

unique platform

for P.A.D.r abutments		
	PADR-INT-CAMTRA-L	PADR-ANA-DIG
for PLAIN B-Space abutments		
	A-INT-CAMTRA-PLAIN-L	A-PLAIN-ANA-330-DIG
for Conico abutments		
	C-INT-CAMTRA-L	ANA2-MD-TS-DIG

Screws and screwdrivers for digital analogs

unique platform

					
VADIG-180	PIN-ANA-DIG	INS-ANA-140	INS-ANA-180	INS-ANA-180-L	INS-ANA-200
screw for digital analog M1.8	pin for digital analog	carrier for digital analog M 1.4 P.A.D.r	carrier for digital analog M1.8	carrier for digital analog M1.8, long	carrier for digital analog M 2.0 Conico

DIGITAL POWER DIGITAL POWER DIGITAL POWER DIGITAL POWER DIGITAL POWER

unique platform post h

repositionable	6.00	 Diagram of the L-MDT-340-6 repositionable drill bit. It shows a yellow drill bit with a diameter of $\varnothing 3.40$ and a length of 6.00. The bit has a double-flute design and a repositionable tip.	L-MDT-340-6
	10.00	 Diagram of the L-MDT-340-10 repositionable drill bit. It shows a yellow drill bit with a diameter of $\varnothing 3.40$ and a length of 10.00. The bit has a double-flute design and a repositionable tip.	L-MDT-340-10
rotating	6.00	 Diagram of the L-MDT-340-6-ROT rotating drill bit. It shows a yellow drill bit with a diameter of $\varnothing 3.40$ and a length of 6.00. The bit has a double-flute design and a rotating tip.	L-MDT-340-6-ROT
	10.00	 Diagram of the L-MDT-340-10-ROT rotating drill bit. It shows a yellow drill bit with a diameter of $\varnothing 3.40$ and a length of 10.00. The bit has a double-flute design and a rotating tip.	L-MDT-340-10-ROT

For these components use the standard screw VM2-180 or the angled hole screw A-VMA-180. Screws are not included.

unique platform cone h 5.00

transgingival h	
	A-TCS-330S-C
1.50 mm	
	A-TCS-330S-C-15
2.50 mm	
	A-TCS-330S-C-25

unique platform

total h 12.00 mm
cementing cone
h 4.00 mm

4.50

ø 4.20

A-CCBAS-330-4


37

T-Connect Unica


unique platform




A-TC-UNI-330




A-TC-UNI-330-1




A-TC-UNI-330-2




A-TC-UNI-330-3




A-TC-UNI-330-ROT



A-TC-UNI-330-1-ROT



A-TC-UNI-330-2-ROT



A-TC-UNI-330-3-ROT

repositionable

rotating

VM2-180 screw included.

ATTENTION: The angled hole screw must always be ordered separately because it is not included in the Interfase package. If you do not have the appropriate screwdriver DSPDCLH-24 or DSPDCLH-32, please order it separately.

Cutting guide and holder for T-Connect Unica



Cutting guide holder S
(small) for T-Connect Unica
left side

TC-UNI-HOLD4-SX

Cutting guide holder S
(small) for T-Connect Unica
right side

TC-UNI-HOLD4-DX


Cutting guide holder L (large)
for T-Connect Unica
right side

TC-UNI-HOLD6-DX

Cutting guide holder L (large)
for T-Connect Unica left side


TC-UNI-HOLD6-SX

4 mm cutting guide for
T-Connect Unica



TC-UNI-DIMA-H4

6 mm cutting guide for
T-Connect Unica



TC-UNI-DIMA-H6

Cutting guide for angled hole for T-Connect Unica



TC-UNI-DIMA-ANG

Please note: STL file available for independently 3D printing of the holders.



See how easy it is to precisely customize the T-Connect Unica with the cutting guides and the corresponding holder.

T-Connect Unica accessories

Intraoral testing O-ring for T-Connect Unica




TC-UNI-ORING

Castable sleeve for
T-Connect, rotating



TC-UNI-CC

Small castable sleeve
for T-Connect Unica,
repositionable



TC-UNI-CC-S

Large castable sleeve
for T-Connect Unica,
repositionable



TC-UNI-CC-L



O-Ring inserter



TC-UNI-ORING-INS


Prama IN - Seating on the neck marginless or with margin

Prama IN T-Connect closing around the implant neck at 0.50 mm

cone height		ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
repositionable	5.00 mm			
		L-TCIN-380-05	L-TCIN-425-05	L-TCIN-500-05
rotating	5.00 mm			
		L-TCIN-380-05-ROT	L-TCIN-425-05-ROT	L-TCIN-500-05-ROT

Only for Regular and Long Neck | For these components use the standard VM2-180 fixation screw or the angled screw hole one A-VMA-180. Screws are not included.

Prama IN T-Connect closing around the implant neck at 1.00 mm, Regular Neck

cone height		ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
repositionable	5.00 mm			
		L-TCIN-380-10	L-TCIN-425-10	L-TCIN-500-10

Only for Regular Neck | For these components use the standard VM2-180 fixation screw or the angled screw hole one A-VMA-180. Screws are not included.

Prama IN T-Connect closing around the implant neck at 1.00 mm, Long Neck

cone height		ø 3.80	ø 4.25	ø 5.00
repositionable	5.00 mm			
		LL-TCIN-380-10	LL-TCIN-425-10	LL-TCIN-500-10

Only for Long Neck | For these components use the standard VM2-180 fixation screw or the angled screw hole one A-VMA-180. Screws are not included.

Cuttable Prama IN T-Connect T-CUT closing around the implant neck at 0.50 mm

cone height		ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
8.00 mm				
		L-TCIN-380-05-8	L-TCIN-425-05-8	L-TCIN-500-05-8
8.00 mm				
		L-TCIN-380-05-8-ROT	L-TCIN-425-05-8-ROT	L-TCIN-500-05-8-ROT

Only for Regular Neck | For these components use the standard VM2-180 fixation screw or the angled screw hole one A-VMA-180. Screws are not included.

Cuttable Prama IN T-Connect T-CUT closing around the implant neck at 1.00 mm, Regular Neck

cone height		ø 3.80 ●	ø 4.25 ●	ø 5.00 ●
8.00 mm				
		L-TCIN-380-10-8	L-TCIN-425-10-8	L-TCIN-500-10-8

Only for Regular Neck | For these components use the standard VM2-180 fixation screw or the angled screw hole one A-VMA-180. Screws are not included.

Cuttable Prama IN T-Connect T-CUT closing around the implant neck at 1.00 mm, Long Neck

cone height		ø 3.80	ø 4.25	ø 5.00
8.00 mm				
		LL-TCIN-380-10-8	LL-TCIN-425-10-8	LL-TCIN-500-10-8

Only for Long Neck | For these components use the standard VM2-180 fixation screw or the angled screw hole one A-VMA-180. Screws are not included.

ATTENTION: The angled hole screw must always be ordered separately because it is not included in the Interfase package. If you do not have the appropriate screwdriver DSPDCLH-24 or DSPDCLH-32, please order it separately.

T-Connect Prama IN accessories



TCIN-HOLD4-SX

TCIN-HOLD4-DX

TCIN-HOLD6-DX

TCIN-HOLD6-SX

TCIN-DIMA-H5

TCIN-DIMA-H6



TCIN-DIMA-ANG

Please note: STL file available for independently 3D printing of the holders.

Ø 3.80 ●

Ø 4.25 ●


Ø 5.00 ●



L-TCS-380S-05



L-TCS-425L-05



L-TCS-500L-05

**Only for Regular and Long Neck. L-VMS screw included.*

unique platform

repositionable



in Titanium

A-CAMT-330-12-SW



in Cobalt Chrome

A-CAMC-330-12-SW

VM2-180 screw included.

Prama IN closing around the neck at 0.50 mm*

Ø 3.80 ●

Ø 4.25 ●

Ø 5.00 ●

repositionable



in Titanium

L-CAMT-38005-12-SW



L-CAMT-42505-12-SW



L-CAMT-50005-12-SW

VM2-180 screw included. | *Only for Regular and Long Neck.

unique platform

repositionable



in Titanium


A-CAMTABU330-CRE

VM2-180 screw included.


Solutions for overdenture

Ball attachments


unique platform



A-AS-330-1



A-AS-330-2



A-AS-330-4

Analog for ball attachments

unique platform



ANAS

Caps for ball attachments

unique platform

polyamide cap



040CRNSNDR8-SM*

steel container for polyamide cap external diameter
ø 4.80 mm



041CAN-SM*

Gr. 5 titanium cap with cap in two parts, retention spring titanium and plastic, plastic mounting ring



CAP-TIT-1

spare plastic ring



AN-CAP-TIT-1

spare steel spring, average hardness



MOL1-CAP-TIT-1

spare steel spring, soft



MOL2-CAP-TIT-1

gold alloy cap, with 3 plastic rings and a tin spacer for the lab



CAP-1

* Pack of 6 pieces.

Overdenture on bars

castable bar, l 5.00 cm, h 3.00 mm, thickness 2.20 mm, ovoid-shaped profile with spacer



BARC-CAV-TIT

divisible bar attachment in titanium for oval bars, h 3.00 mm, thickness 2.20 mm



CAV-TIT

castable bar, l 5.00 cm, ø 2.20 mm



BARC

bar attachment in gold alloy for round bars ø 2.20 mm



CAV-375

Locator abutment

transgingival h

unique platform

1.00 mm



1670

2.00 mm



1671

3.00 mm



1672

4.00 mm

-

5.00 mm



1674

6.00 mm

-

Accessories for Locator abutment

unique platform

4 aluminium analogs



8530

4 aluminium transfers,
4 retainers cod. 8515 included



8505

4 polyethylene parallelism pins



8517

steel plate AISI 316 L



9530

complete Locator Core Tool



8393

steel push rod for
inserting the retainers



8397

steel driver for abutment
screwing/unscrewing



8390

retention jacket for the driver, to
transfer the Locator abutments
into the oral cavity



8394

Retainers and spacer rings for Locator abutments

unique platform

low retention capacity for impression taking,
black, in polyethylene,
4 pcs packaging



8515

retention 5 lb = 2268 g, in nylon,
disparallelisms up to 10° per side,
4 pcs packaging



8524

retention 3 lb = 1361 g, in nylon,
disparallelisms up to 10°
per side, 4 pcs packaging



8527

retention 1,5 lb = 680 g, in nylon,
disparallelisms up to 10°
per side, 4 pcs packaging



8529

retention 4 lb = 1814 g, in nylon,
disparallelisms up to 20°
per side, 4 pcs packaging



8547

retention 2 lb = 907 g, in nylon,
disparallelisms up to 20°
per side, 4 pcs packaging



8915

retention 1 lb = 450 g, in nylon,
disparallelisms up to 20°
per side, 4 pcs packaging



8548

spacer rings in silicon rubber,
20 pcs packaging



8514

kit: 2 Gr. 5 titanium caps, 2 spacer rings,
2 black retainers and 2 nylon retainers for
each of the 3 different retention capacities



8519-2

kit: 2 Gr. 5 titanium caps, 2 spacer rings,
2 black retainers and 2 nylon retainers for
each of the 3 different retention capacities
meant for disparallelisms



8540-2

kit: 2 Gr. 5 steel caps, 2 spacer rings,
2 black retainers and 2 nylon retainers for
each of the 3 different retention capacities



8550-2

Locator R-Tx abutment*

transgingival h	unique platform
1.00 mm	 31110-01
2.00 mm	 31110-02
3.00 mm	 31110-03
4.00 mm	 31110-04
5.00 mm	 31110-05
6.00 mm	 31110-06

* All-in-One packaging that includes an abutment, a spacer ring, a cap, and an insert for each retention (gray, blue, pink, white).

Caps and retainer for Locator R-Tx abutment

	unique platform
pack of 4 caps and 4 black retainers for impression taking	 30013
grey, without retention, pack of 4 pcs	 30001
blue, low retention, pack of 4 pcs	 30002
pink, medium retention, pack of 4 pcs	 30003
white, high retention, pack of 4 pcs	 30004
black, for the impression taking, pack of 4 pcs	 30012

Accessories for Locator R-Tx abutment

	unique platform
aluminium analogs for abutments $\varnothing$ 3.35 mm pack of 4 pcs	 30014
aluminium analogs for abutments $\varnothing$ 4.00, pack of 4 pcs	 30015
aluminium analogs for abutments $\varnothing$ 5.00, pack of 4 pcs	 30016
transfers with low retention, pack of 4 pcs	 30017
spacer rings for relining phase, pack of 4 pcs	 30018
Locator core tool for Locator R-Tx	 30021
RA screwdriver for Locator R-Tx, hex 1.25 mm	 08008
digital screwdriver for Locator R-Tx, hex 1.25 mm	 08007

Clinical and laboratory prosthetic screws

Standard fixation screw		
	VM2-180	
	VM2-180-10	10 pcs packaging
Conical support fixation screw		
	L-VMS-180	
	L-VMS-180-10	10 pcs packaging
Full Head Fixation screw		
	L-VM-180	
	L-VM-180-10	confezione da 10 unità
Fixation screw for angled screw hole		
	A-VMA-180	
Screws for Pick-up transfer		
	VTRA2-180-15	h 15.00
	VTRA2-180-20	h 20.00
	VTRA2-180-MAN	for manual screwing
Screw for closed-tray transfers		
	VTRA2-180-10	h 10.00
Screws for P.A.D.r abutments		
	PADR-VM-180	for P.A.D.r abutment
	PADR-VM-180-10	for P.A.D.r abutment 10 pcs packaging
	PAD-VTRAL-140	for Pick-up transfer, long
	PAD-VTRA-140	for Pick-up transfer, standard
	PAD-VTRAL-140-MAN	for Pick-up transfer, manual screwing
	PAD-VP-140	for prosthetic components on P.A.D.r abutments
	PAD-VP-140-10	for prosthetic components on P.A.D.r abutments 10 pcs packaging
	PAD-VPA-140	fixation screw for PADR-TCA-M to use with the dedicated screwdriver DSPDCLH-24 and DSPDCLH-32
	PADR-VTRAL-140-M	carrier for P.A.D.r

Screws for CAD-CAM prosthesis		
	A-PLAIN-VP200	for PLAIN abutment
	A-PLAIN-VP200-10	for PLAIN abutment, 10 pcs packaging
	A-PLAIN-VTRA200	for PLAIN abutment Pick-up transfers
	ABU-VMA-200	for angled screw hole
Screws for CAD-CAM prosthesis		
	A-CAMTVABU180	for zirconium oxide single crowns and ECHO overstructures
	PAD-VCAM140	for zirconium oxide overstructure on abutment P.A.D.r
	CAMPRON205-10	spare ring washers for the head of the screw
	VADIG-180	locking screw for digital analogs

Recommended torque for tightening prosthetic components	
8-10 Ncm	surgical cover screws, healing abutments, transfer screws
20-25 Ncm	prosthetic screws
25-30 Ncm	direct screwing on implants prosthetic components
8-10 Ncm	transfer screws

Warning

Fixation screws are always included in the packaging of each prosthetic component, unless otherwise indicated. Please use laboratory screws for laboratory phases and keep the screw included with the prosthetic component for the final fastening in the patient's mouth.

Bibliography

- a. Canullo L., Tallarico M., Pradíes G., Marinotti F., Loi I., Cocchetto R.,
Soft and hard tissue response to an implant with a convergent collar in the esthetic area: preliminary report at 18 months,
Int J Esthet Dent 2017; 12:2-19:
- b. Canullo L, Menini M, Covani U, Pesce P,
Clinical outcomes of using a prosthetic protocol to rehabilitate tissue-level implants with a convergent collar in the esthetic zone: a 3-year prospective study,
J Prosthet Dent. 2019 Jun 18. pii: S0022-3913(19)30077-0. doi: 10.1016/j.prosdent.2018.12.022
- c. Valente N.A., Wu M., Toti P., Derchi G., Barone A,
Concave/convergent versus parallel/divergent implant transmucosal profiles impact on hard and soft peri-implant tissues: a systematic review with meta-analysis,
Int J Prosthodont. Sep/Oct 2020;33(5):553-564. doi: 10.11607/ijp.6726
- d. Cabanes Gumbau G, Pascual-Moscardò A, Penarrocha-Oltra D, Garcia-Mira B, Aizcorbe-Vicente J, Penarrocha-Diogo M,
Volumetric variation of peri-implant soft tissues in convergent collar implants and crowns using the biologically oriented preparation technique,
Med Oral Patol Oral Cir Bucal. 2019 Sep 1;24(5):e643-e651. doi: 10.4317/medoral.22946
- e. Agustín-Panadero R, Martínez-Martínez N, L. Fernandez-Estevan, J. Faus-Lopez, F. Sola-Ruiz,
Influence of transmucosal Area Morphology on Peri-implant bone loss in tissue-level implants,
The International Journal of Oral & Maxillofacial Implants, 2019; 34: 947-952. doi: 10.11607/jomi.7329
- f. Marchio V, Derchi G, Cinquini C, Miceli M, Gabriele M, Alfonsi F, Barone A,
Tissue level implants in healthy versus medically compromised patients: a cohort comparative study,
Minerva Stomatologica 2020 mese;69(0):000–000, Doi: 10.23736/S0026-4970.20.04359-9
- g. Prati C, Zamparini F, Canullo L, Pirani C, Botticelli D, Gandolfi MG,
Factors affecting soft and hard tissues around two-piece transmucosal implants: a 3 year prospective cohort study,
J Oral Maxillofac Implants 2020;35, doi: 10.11607/jomi.7778
- h. Agustín-Panadero, R.; Bermúdez-Mulet, I.; Fernández-Estevan, L.; Fernanda Solá-Ruiz, M.; Marco-Pitarch, R.; García-Selva, M.; Zubizarreta-Macho, Á.; León-Martínez, R.,
Peri-Implant Behavior of Tissue Level Dental Implants with a Convergent Neck.,
Int. J. Environ. Res. Public Health 2021, 18, 5232, Doi: 10.3390, ijerph18105232
- i. Palombo D., Rahmati M., Vignoletti F., Sanz-Esporrin J., Haugen H. J., Sanz M.,
Hard and soft tissue healing around implants with a modified implant neck configuration: An experimental in vivo preclinical investigation,
Clin Oral Impl Res. 2021;00:1–15.
- j. Morón-Conejo B, Sanz-Sánchez I, Salido MP, Martínez-Rus F, Pradíes G,
The effect of a convergent transmucosal neck on soft tissues and radiographic outcomes: a 1-year follow-up randomized controlled trial,
Clin Oral Investig 2023 Feb 7. doi: 10.1007/s00784-023-04892-9.
- k. Mandillo-Alonso V, Cascos Sanchez R, Antonaya-Martin JL, Laguna-Martos M,
Evaluation of peri-implant soft and hard tissues behavior in screw-retained crowns by the biologically oriented preparation technique: Ambispective longitudinal analytical study,
J Clin Exp Dent. 2022 Jan; 14(1): e64–e71, doi: 10.4317/jced.58924
- l. Mandillo-Alonso V, Cascos Sanchez R, Antonaya-Martin JL, Laguna-Martos M,
Soft tissue thickness evaluation in screw-retained crowns by the biologically oriented preparation technique,
J Clin Exp Dent 2021 Dec 1;13(12): e1209-e1215, doi: 10.4317/jced.58952



rev. 05-25



**Authorized Distributor
Sweden & Martina Inc.**

78 John Miller Way
Unit 1021 - Kearny, New Jersey 07032
Toll free 1-844-8MARTINA (1-844-862-7846)
info.us@sweden-martina.com

www.sweden-martinainc.com

The implants, standard prosthetic components and surgical instruments contained in this issue are medical devices and are manufactured by Sweden & Martina S.p.A. except for:

- the Locator Abutment and related components which are medical devices manufactured by Zest Anchors, Inc., 2875 Loker Avenue East, Carlsbad, CA 92010, USA. Locator is a trademark by Zest Anchors, Inc.
- The TWL torque wrench with control lever and its adapter are medical devices manufactured by Elos Medtech Pinol A/S, Engvej 33, 3330 Gorlose, Denmark.

They conform to the ISO 9001 and ISO 13485 standards and are certified with the CE in compliance with Regulation (EU) n. 2017/745 for medical devices.
We have met the good manufacturing standards (GMP) set forth by many countries worldwide, including the United States FDA.



Some products may not be regulatory/released for sale in all markets.

All trademarks herein are the property of Sweden & Martina S.p.A. unless otherwise indicated.

This material is intended for laboratories and clinicians and is not intended for patient distribution.

This material is not to be redistributed, duplicated, or disclosed without the express written consent of Sweden & Martina S.p.A.

For additional product information, including indications, contraindications, warnings, precautions, and potential adverse effects, see Sweden & Martina S.p.A. website.

The contents are updated at the time of publication. Check with the company for any subsequent updates.